Occlusive Design	Patented, woven wire design prevents the migration of fragments as small as 1 mm during intracorporeal lithotripsy. ¹
Stone Release	Flexible “umbrella design” readily releases stones too large for extraction. ²
145 cm Length	145 cm length allows for expeditious scope exchange after placement under direct vision.
Nitinol Alloy	Resilient memory characteristics retain preformed shape when deployed.


3

After reintroduction of the endoscope through the access sheath, an intracorporeal lithotripter may now be utilized to fragment the stone. Maintain a safe distance between all intracorporeal lithotripters and the NTrap during stone fragmentation. Care should be taken when utilizing Holmium Lasers to avoid firing the laser directly on the NTrap.


4

During intracorporeal lithotripsy, the NTrap prevents migration of stones and stone fragments as small as 1 mm.¹

5

After successful fragmentation of the stone, the NTrap is utilized to entrap remaining stone fragments for extraction. The highly flexible nitinol wires are designed to readily and safely disengage from any stone fragment that is too large for safe extraction.² Again, a ureteral access sheath may be used to provide a protected conduit for instrument and stone removal.

NTrap Nitinol Stone Extractor

NTrap is an endoscopic entrapment and extraction device for calculi and other foreign bodies in the urinary tract, effectively minimizing stone migration during laser, ultrasonic, electrohydraulic or pneumatic lithotripsy.

Global Product Number	Order Number	Sheath French Size	Sheath Tip French Size	Length cm	Configuration
G32724	NTP-028145	2.8	3.0	145	7 mm umbrella design

www.cookmedical.com

Cook Europe Shared Service Centre

Country		Telephone	Email
Austria		(+43) 1-795-67-121	oeorders@cook.ie
Belgium	French	(+32) 27-00-16-33	beorders@cook.ie
Belgium	Flemish	(+32) 27-00-16-33	nlorders@cook.ie
Denmark		(+45) 38487607	daorders@cook.ie
Distributors		(+353) 61-239240	sscdistributors@cook.ie
France		(+33) 1-71-23-02-69	frorders@cook.ie
Germany		(+49) 6950-072804	deorders@cook.ie
Italy		(+39) 02-69-68-28-53	itorders@cook.ie
Netherlands		(+31) 202-013-367	nlorders@cook.ie
Norway		(+47) 23-16-2968	noorders@cook.ie
Spain		(+34) 91-2702691	esorders@cook.ie
Switzerland	French	(+41) 44-8-009-609	frorders@cook.ie
Switzerland	Italian	(+41) 44-8-009-609	itorders@cook.ie
Switzerland	German	(+41) 44-8-009-609	deorders@cook.ie
Sweden		(+46) 85-876-94-68	seorders@cook.ie
UK		(+44) 20-7365-4183	ukorders@cook.ie

AORTIC INTERVENTION

CARDIOLOGY

CRITICAL CARE

ENDOSCOPY

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH

NTrap[®]

NITINOL STONE EXTRACTOR

NTrap[®]

NITINOL STONE EXTRACTOR


NTrap continues the groundbreaking series of nitinol designs from Cook Medical, the world leader in nitinol basket technologies.

When deployed in the ureter, the NTrap's intricately woven, multi-wire design provides an occlusive umbrella, preventing retrograde stone migration during intracorporeal lithotripsy.

1

At deployment, the NTrap should be positioned under direct vision through the endoscope. Once the NTrap is safely beyond the ureteral stone, it is deployed to occlude the ureter.

2

Leaving the NTrap deployed above the stone, the endoscope may now be easily removed to free up the working channel and optimize irrigant flow. Utilization of a previously placed ureteral access sheath* will ease this process and protect the ureter during instrument exchanges.

* Flexor[®] Ureteral Access Sheath, Cook Medical

References:

1. Holley PG, Sharma SK, Perry KT, et al. Assessment of novel ureteral occlusion device and comparison with stone cone in prevention of stone migration during lithotripsy. *J Endourol.* 2005;19(2):200-3.
2. Sharma SK, Ouwenga M, Holley P, et al. Load-release points of two novel ureteral stone-trapping devices. *J Endourol.* 2005;19(7):894-7.

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