H-30™

HOLMIUM LASER SYSTEM

CE
0088

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Table of Contents

Chapter 1: Introduction

Introduction to the H-30™	1-4
About Cook Medical	1-4
Holmium Lasers	1-5
Laser Lithotripsy	1-5
Device Description	1-5 - 1-7
System Specifications	1-8
Labels and Symbols	1-9

Chapter 2: Installation

Installation	1-11
Transportation	1-11
Unpacking the Laser System	1-11 - 1-12
Contents	1-12
Ambient Operating Conditions	1-12
Filling the Cooling System with Water	1-12 - 1-13
Initial System Power On	1-13 - 1-14
Verifying Alignment and Calibration	1-14

Chapter 3: Safety

General Warnings and Caution	1-16
General Safety Information	1-16
Staff Training	1-16
Working Area	1-16
Repair and Modifications of the Device	1-16
Eye and Skin Exposure	1-17
Hazards	1-17 - 1-18
Interference with Other Devices	1-18

Chapter 4: Indications and Instructions for Use

Indications for Use	1-20
Surgical Applications	1-20
Warnings and Risks	1-20 - 1-21
Parameter Determination	1-21
Setting the Pulse Width	1-21 - 1-22

Adjusting the Energy per Pulse	1-22
Surgical Complications and Risks	1-22 - 1-23
Contraindications for Laser Surgery	1-23

Chapter 5: Operation

Laser Screens	1-25
Main Screen	1-25 - 1-26
Options Screen	1-27 - 1-28
Laser Parameters	1-28 - 1-29
Turning the System Off	1-29
Protection Against Unauthorized Use	1-29
Error Messages	1-29 - 1-30

Chapter 6: Maintenance

General Maintenance	1-32
Cooling Water Level	1-32
Power Cord	1-32
Labels	1-32
Replacing the Blast Shield	1-33 - 1-34
Troubleshooting	1-35
Error Messages	1-35 - 1-36

Chapter 7: Warranty Information

Warranty Statement	1-38
Limits and Exclusions	1-38
Warranty Period	1-39
Warranty Procedure	1-39

Chapter 1

Introduction

[Return to Table of Contents](#)

Introduction to the H-30™

The Cook Medical H-30™ Holmium Laser System and holmium laser fibers are effective in the treatment of soft tissue disease and urinary stones. The H-30 has indications in a variety of medical specialties including urology.

The following are urological indications for use: urethral strictures, bladder neck incisions, ablation, and resection of bladder tumors, urethral tumors and ureteral tumors.

Multiple uses for lithotripsy and percutaneous urinary lithotripsy are also indicated: endoscopic fragmentation of urethral, ureteral, bladder and renal calculi, endoscopic fragmentation of kidney calculi, and the treatment of distal impacted fragments of steinstrasse when a wire guide cannot be passed.

About Cook Medical

Choosing Cook Medical means choosing quality and innovation without compromise. To attain the best results with the H-30 and avoid the risk of dangerous faults, please thoroughly read this user manual before using the laser system.

Contact your sales representative if you have any questions regarding the operation or performance of the H-30 Holmium Laser System and holmium laser fibers.

Cook Medical

750 Daniels Way, P.O. Box 489
Bloomington, IN 47402-0489 USA
800.457.4500 ext. 102146
customerrelationsNA@cookmedical.com

www.cookmedical.com

Holmium Lasers

Holmium (Ho:YAG) lasers are versatile, covering a wide range of indications. Holmium laser energy is largely absorbed by water, making it ideal for tissue ablation with minimal lateral thermal damage. In the case of laser lithotripsy, laser energy vaporizes water in the stone, causing it to break into smaller particles that pass easily through the urinary tract.

Laser Lithotripsy

Holmium lasers can treat any type of calculi, regardless of color or composition, from the bladder to the kidney. A major advantage of a holmium laser is that its laser fibers can be placed through small, flexible ureteroscopes with minimal loss of deflection or irrigation. Because of the versatility and safety of holmium laser energy, the indications for treatment have expanded to include stones larger than 2 cm.

Device Description

Name	Laser Source
H-30 Holmium Laser System	Ho:YAG

Front of Laser

Figure 1

Top Rear of Laser


Figure 2

Bottom Rear of Laser


Figure 3

Device Connections

The H-30 contains one laser source, emitting at 2.1 μm . The laser is delivered through an optical fiber (holmium laser fiber) connected to the Laser Aperature. (See figure 1.)

Key Switch

The key switch turns the device on. There are two positions on the switch: 'I' (on) and 'O' (off). To power on the device, insert the key and turn it clockwise to the on position. To power off, turn the key to the off position and remove it. (See figure 2.)

Emergency Stop Button

The emergency stop button allows the immediate shutdown of the device. In case of an emergency, push the emergency stop button. To reset the emergency stop button, twist and pull it out. (See figure 1.)

Control Panel


Figure 4


System Specifications

General Specifications	
Electrical Requirements	230 V ~ ±10%/50-60 Hz/10 A (30 W output) 115 V ~ ±10%/50-60 Hz/20 A (30 W output) 115 V ~ ±10%/50-60 Hz/15 A (20 W output)
Dimensions	20 in (W) x 23 in (D) x 51 in (H) 50 x 59 x 129 cm
Weight	230 lbs 104 kg
Operation Temperature	50°F - 77°F, 10°C - 25°C
Storing Temperature	41°F - 104°F, 5°C - 40°C
Relative Humidity	From 30% to 85% (no condensing)
Cooling	Air cooled
Electrical Protection Class	I
Type	B
Degree of Protection	IP20 (IEC529)
Laser Specifications	
Laser Type	Flashlamp pumped CTH:YAG
Wavelength	2.1 μm
Power to Tissue	30 W
Energy/Pulse	0.5 - 3.5 J
Application Mode	Pulsed
Pulse Width	Short and Long
Repetition Rate	5 - 20 Hz
Delivery Device	Wide range of flexible silica fibers
Aiming Beam	Diode laser, green, 532 nm (adjustable < 3 mW)

Labels and Symbols


WARNING—Visible and invisible laser radiation.

Laser Stop

Laser Aperture

Optical Fiber Applicator

Foot switch
laser activation device

Chapter 2

Installation

[Return to Table of Contents](#)

Installation

The H-30 laser system should be installed by trained personnel upon initial receipt of the laser system from Cook Medical. Any damage that occurs after transportation must be corrected prior to use.

WARNING! Do not attempt to use the H-30 prior to trained personnel performing the installation procedure.

Transportation

Transportation of the laser system must be performed in compliance with all local safety laws.

The laser system must be packaged to protect it from damage, such as in the original shipping crate. If the original shipping crate is not used for transportation, the system must be secured to prevent damage. The handles can be used to tie the system to a secure surface. The crate should only be moved using proper lifting and transporting means, such as a pallet jack.

The laser system should not be subjected to excessive force, shock or vibration. If excessive shock or vibration is suspected to have occurred, the laser alignment should be inspected by a trained technician.

If any water remains in the cooling system, the laser must not be subjected to freezing temperatures. Freezing water within the laser cavity can damage optical components such as the laser cavity, rod and flash lamp or mechanical components such as the heat exchanger or water pump. If there is a possibility of freezing temperatures, the laser's cooling system should be drained of all water by attaching fill/drain tubes to both the Drain and Air Purge ports on the rear panel of the laser system. (See figures 2 and 3.)

Storage humidity should not exceed 95%. If the laser is stored in a high humidity location, it must be allowed sufficient time to reach equilibrium when introduced to a cool operating environment. Failure to do this may result in damage to the optics from condensation.

Unpacking the Laser System

Your H-30 laser system is shipped from the factory in a wooden crate to protect it from damage during shipping. Before removing the laser system, inspect the crate for damage. If there is evidence of damage, save all packing materials and notify the shipping company by filing an insurance claim report.

A Cook sales or service representative will help you with the installation of the laser system. If you have any questions, please contact Customer Relations. After the laser system has been unpacked, please store the shipping crate and packing materials in a safe location in case the laser system has to be transported in the future.

The crate can easily be disassembled and laid flat to be stored more easily.

Contents

The H-30 laser system is shipped with the following components and accessories. Please take inventory of the items listed below and contact Customer Service if any items are missing or damaged.

- Laser system
- Two keys to key switch
- Foot switch
- Remote interlock connector
- Container of distilled water
- Plastic water bottle for filling the cooling system
- Two fill/drain/air purge tubes

Ambient Operating Conditions

The laser system should be used in a well-ventilated, air-conditioned operating room. The procedure room temperature should be kept between 18-24°C (65-75°F) for optimum system performance. In general, holmium lasers are very sensitive to temperature. When the laser system reaches above 38°C (100°F), the safety circuits will prevent the laser system from operating. The information bar on the touchscreen will read "Over temperature." The laser system may overheat for the following reasons:

- Improper ventilation
- The ambient operating room temperature is above 24°C (74°F)
- The laser system has been running at high power for an extended period of time (more than 30 minutes)
- The input line voltage is lower than normal, causing the cooling fans or water pump to under-perform

Filling the Cooling System with Water

1. Connect a fill/drain tube to the laser system in the port labeled Air Purge.
2. Fill the supplied water bottle with distilled water and connect its output tube to the water fill port on the laser system labeled Water Fill. (See figure 2.)

NOTE: Use distilled water only.

3. Direct the output of the Air Purge tube into the distilled water container used to fill the water bottle.
4. Hold the water bottle up and compress the output valve to allow water to flow into the laser system.
5. Once some water has started to flow into the laser system, turn the laser system on and continue filling the cooling system.
6. Turn off the laser system and refill the water bottle whenever necessary.
7. Repeat the process until only water is forced out of the Air Purge port.

NOTE: The laser system has built-in safety circuitry to detect when there is too little water flow. If this circuitry engages while it fills with water, the information bar will display "Low coolant flow." Simply turn the laser system off and back on to continue.

8. Disconnect both the Water Fill and Air Purge tube and let the laser system run for a few minutes to force any air in the system to the top.
9. Inspect the water level indicator on the rear panel. Repeat procedure until the red float is at the top of the level indicator. (See figure 2.)
10. Store tubes and water bottle in a safe place for future use.

Initial System Power On

Review the previous sections in this chapter and verify that all site installation requirements have been met. The laser system should only be operated in an area of restricted access. Post a laser warning sign on each access door.

1. Connect the foot switch and remote interlock connector to the laser system. (See figure 3.) Find the remote interlock connector supplied separately in the bag that comes with the unit.
2. Connect the laser device to a suitable AC electrical outlet.
3. Ensure the emergency stop button isn't engaged. (See figure 1.)
4. Turn the MAINS switch on. (See figure 3.)
5. Turn the key switch on. (See figure 2.)
6. Lift the laser aperture cover and connect a 550 μm holmium laser fiber to the laser. (See figure 1.)

NOTE: Do not over-tighten the laser fiber. The connector only needs to be finger tight for proper operation.

7. After a few seconds, the information bar on the control panel screen will display "STANDBY." (See figure 4.)
8. Press the Ready button. (See figure 4.) The Ready progress bar will take about two seconds to fill.

9. Once the Ready bar is fully lit, the information bar will display "READY to lase." This indicates that the laser system is operating normally. (See figure 4.)
10. Follow the instructions in the next section to verify alignment and calibration.

Verifying Alignment and Calibration

The laser system has been designed with many safeguards to prevent optical misalignment due to shipping and normal operation. However, it is important to verify the alignment and calibration of the laser system to ensure peak system performance. The following test requires an energy or power meter calibrated to 2100 nm.

The following test is intended for laser systems capable of delivering the full 30 W output power.

1. Set up laser system for lasing.
2. Verify that the 550 μm test fiber has clean surfaces on both the proximal connector and distal ends.
3. Insert the 550 μm test fiber into the laser aperture.

WARNING! Wear appropriate safety eyewear. While the laser is in operation, all personnel in the operating room must wear appropriate laser eyewear.

4. Set up laser power probe in front of and parallel to the distal face of the test fiber, as shown. Be careful not to place the fiber tip too close to the probe, as the energy density may damage the probe's sensitive surface.
5. Set the power parameters to Short pulse width mode, 10 Hz, 0.5 J.
6. Press the Ready button and wait until the information bar displays "READY to lase."
7. Depress the foot switch for 5-10 seconds until you attain a stable reading on the power meter.
8. If the reading is not within 10% of the displayed settings, the laser system is either misaligned or out of tune. Please contact your sales representative.
9. Repeat steps 7 through 9 for the following settings as well:
 - a. Short pulse width mode, 10 Hz, 2.0 J.
 - b. Short pulse width mode, 10 Hz, 3.0 J. (NOTE: Only for max output of 30 W.)


Figure 5

Chapter 3

Safety

[Return to Table of Contents](#)

General Warnings and Caution

The H-30 is a 30 W holmium laser system for surgical procedures. The device is manufactured following the 93/42/EEC Annex II standard and it is identified as follows:

Device	30 W Holmium Laser System
Name	H-30 Holmium Laser System
Category of the Device	Laser for medical use
Manufactured By	Cook Urological Inc. 1100 West Morgan Street Spencer, IN 47460-9426 USA

General Safety Information

- All persons operating this equipment must understand the instructions for operation and the instructions for safety specified in this manual.
- Only authorized individuals with appropriate medical laser training and knowledge should operate the laser system.
- Only authorized service personnel should access the internal electrical components.
- Access to the user manual should be available in the working area of the laser.
- All the warning labels must be kept in good condition.

DO NOT dispose of this laser equipment in any location other than designated locations.

Staff Training

The use of the laser is restricted only to medical staff with experience using lasers.

Working Area

This H-30 is a Class 4 laser and must be used in a specific working area defined and delimited following local and international standards (IEC/EN 60825-1: 2003-02).

Repair and Modifications of the Device

Only authorized service personnel should execute repairs and maintenance. Constructive modifications are not permitted.

Note: Unauthorized repair or maintenance, and/or modifications to the laser system, will void the warranty.

Eye and Skin Exposure

The working laser beam of the H-30 is invisible and, if exposed to the eyes, may cause sight loss. The skin can resist a significantly higher amount of laser energy, but the laser beam can seriously burn the skin.

A danger exists in any case of:

- Direct laser radiation
- Reflected laser radiation
- Diffused laser radiation

WARNING! All the personnel present in the laser working area must wear appropriate laser safety eyewear.

Use protective eyewear with either of the following specifications:

- I 2100 L2 (EUR) @ 2100NM
- OD 4 (US) @ 2100NM

Always check the integrity and condition of safety eyewear.

Before wearing protective eyewear, be sure that the protective glass and frame are in good condition.

The worst-case nominal ocular hazard distance (NOHD) for the H-30 laser system emitting laser energy from a fiber is 0.5 m. During service operations, there is access to laser energy emitting directly from the laser source, which raises the NOHD to 15 m.

Hazards

Electrical Hazards

- Only a qualified service technician may remove the enclosure covers. Removing the enclosures exposes personnel to potential electrical hazards.
- Do not touch areas marked HIGH VOLTAGE. These and surrounding components of the laser system contain charge voltages of up to 700 VDC.
- Do not operate the laser system if the operating room floor is flooded.
- Service technicians, do not power on the laser system unless all high voltage areas are dry and not in danger of flooding.
- Inspect power cord before use. Do not use a damaged power cord or a power cord with a damaged or loose power plug.

Other Hazards

- **RISK OF FIRE**—Do not point the distal end of the fiber toward flammable materials such as surgical drapes.
- **RISK OF EXPLOSION**—Do not operate the laser system in the presence of flammable substances, including flammable anesthetics.
- **RISK OF UNINTENDED MOVEMENT**—The rear wheels have a directional lock to allow the laser to be easily rolled from room to room when the lock is engaged and easily positioned within an OR when the lock is disengaged. The front wheels on the laser can be locked in place with brakes to keep the laser from shifting during operation.

Interference with Other Devices

Electromagnetic fields from other irradiants in the area surrounding the laser system can affect the laser's performance. The laser's electromagnetic field can also have an affect on other electrical devices.

Do not use cellular phones in the immediate vicinity of the laser system.

Chapter 4

Indications & Instructions for Use

Indications for Use

This section discusses the use of the H-30 Holmium Laser System in clinical applications. The information is provided by specialty and includes procedural recommendations, along with specific indications and contraindications. The information provided in this section is not all-inclusive and is not intended to replace physician training or experience.

Only physicians and staff who have been appropriately trained and who are thoroughly familiar with the instructions and safety precautions provided in this manual should operate the H-30. A review of prior training is strongly encouraged and recommended.

Surgical Applications

The H-30 Laser System and holmium laser fibers are intended for use in surgical procedures using open, laparoscopic and endoscopic incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft tissue.

Urology

Open and endoscopic surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) including:

- Urethral strictures
- Bladder neck incisions (BNI)
- Ablation and resection of bladder tumors, urethral tumors and ureteral tumors
- Endoscopic fragmentation of urethral, ureteral, bladder and renal calculi
- Treatment of distal impacted fragments of steinstrasse when wire guide cannot be passed

Warnings and Risks

General warnings and precautions for the H-30 Holmium Laser System

The physician should become fully acquainted with the unique surgical and therapeutic effects produced with the 2.1 μm wavelength before using the H-30 clinically. These effects include coagulation, depth of penetration and cutting intensity.

Unless otherwise stated in the specific application section, the physician should begin at the lowest power and use shorter durations. The surgeon should note the surgical effect and adjust the settings until the desired result is obtained.

Warnings and precautions for the safe use of holmium laser fibers

- Changes in the tissue texture and color are the best indicators of the laser's effect. The recommended power settings are less important than the observed tissue effect. Specific pulse width depends on the tissue and is left to the surgeon's preference and best medical judgment.
- Use the lowest possible power settings required to achieve the desired tissue effect for treatment.
- Higher energy may be necessary to achieve the desired tissue effects if cooled irrigant is utilized. Higher power settings may increase chances of damage to the laser fiber.
- The use of mechanical pressure on the laser fiber does not increase its cutting or vaporization effects but may induce bleeding, thermal damage and fiber destruction.
- The risk of combustion, perforation, and laser induced hemorrhage, any of which can cause patient death, must be fully explained to the patient prior to surgery.
- Endoscopic treatment may result in adverse reactions such as fever, chills, sepsis, edema, and hemorrhage. In extreme cases, death may occur due to procedural complications, concurrent illness, or the application of laser radiation.

Parameter Determination

Choosing the ideal laser fiber

Cook Medical's holmium laser fibers for use with the H-30 come in 150, 200, 273, 365, 550 and 940 μm core diameters. Larger core fibers will allow the use of higher power but are less flexible. Higher pulse energy and larger fibers result in larger ablation volumes.¹ Large core fibers will also be consumed less during procedures because of the greater thermal mass and ability to dissipate parasitic thermal energy faster. Users are encouraged to select the largest core size feasible considering the flexibility needed to reach the target and the size of the endoscope's working channel.

¹ Kang HW, Lee H, Teichman JMH, et al. Dependence of calculus retropulsion on pulse duration during Ho: YAG laser lithotripsy. *Lasers Surg Med.* 2006;38(8):762-772.

Setting the Pulse Width *(via the toggle button on the main screen)*

Lithotripsy

- Short pulse width will cause larger fragmentation particulates (bust stone) but will also cause greater retropulsion.
- Long pulse width will cause smaller fragmentation particulates (dust stone) and will cause less retropulsion.

Soft Tissue

- Use the Short pulse width if you want to minimize coagulation or create a large, thick cut/vaporization area.
- Use the Long pulse width if you want to increase coagulation and decrease bleeding.

Notes

- Most soft tissue cases will benefit from using the Long pulse width mode.
- If conditions have changed, such as a piece of an impacted stone breaks off or the soft tissue is bleeding more than expected, the pulse width mode can be easily switched with the tap of a toggle button.

Adjusting the Energy per Pulse

Adjust the energy per pulse to match the desired clinical effect.

Lithotripsy

- Lower pulse energy will cause the particulate size to be smaller (dusting the stone) and cause less retropulsion.
- Higher energy will create larger fragments (busting or drilling the stone) but will also cause more retropulsion.

Surgical Complications and Risks

The same complications and risks that exist for conventional or traditional surgery exist for laser surgery such as but not limited to:

- **Pain**—Short lived pain may occur immediately following endoscopic laser surgery and may persist for as long as 48 hours.
- **Fever and Leukocytosis**—Immediately following endoscopic laser surgery, the patient may experience fever and leukocytosis, which are commonly associated with tissue degradation. These conditions are generally resolved without treatment. Tissue cultures are indicated to exclude the possibility of infection.
- **Bleeding**—Patients may experience bleeding at the surgical site or from erosion of a tumor during or after laser surgery. Post treatment blood studies, such as hematocrit levels may be necessary.
- **Sepsis**—Lase ablated tissue may become infected after surgery. If a question of sepsis exists, a culture should be taken and other appropriate evaluations made.
- **Perforation**—Perforation can occur as a result of endoscopic or cystoscopic procedures. To clinically diagnose perforation, patients must be monitored postoperatively through physical signs, hematocrit and radiography.

The following complications can cause patient death:

- **Non-thermal**–Perforation, aspiration, induced hemorrhage, allergic reaction to medication, hypertension, arrhythmia, pain and infection.
- **Thermal Accute**–Induced hemorrhage, ulceration, perforation, edema, pain, fever, leukocytosis and chills.
- **Thermal Chronic**–Delay in healing, perforation, delayed hemorrhage and sepsis.

Contraindications for Laser Surgery

The H-30 should only be used in conditions where its use is appropriate and its efficacy has been proven. A qualified physician should perform clinical applications.

The use of the laser is contraindicated for patients:

- Whose general medical condition contraindicates surgical intervention
- Where appropriate anesthesia is contraindicated by patient history
- Where tissue (especially tumors) is calcified
- Requiring hemostasis of vessels over approximately two millimeters in diameter
- Where laser therapy is not considered the treatment of choice
- Where any irrigant other than physiological saline is used

Chapter 5


Operation

Laser Screens

The touchscreen display contains two functional screens, the main screen and the options screen. Touch the tabs at the top right-hand corner of the display to toggle between these screens.

Main Screen

The main screen allows you to adjust functions for the main power settings.

Main Screen Items (See figure 6)

Main Screen Tab

- Displays power parameters and controls for adjustment
- Touch to access the main screen

Energy Bar

- Displays set pulse energy and total delivered energy (TOTAL) in Joules
- Touch to highlight orange and adjust energy per pulse with the adjustment controls

Frequency Bar

- Displays pulse frequency in Hertz
- Touch to highlight orange and adjust pulse frequency with the adjustment controls

Preset Buttons

- Touch one of the numbered presets to jump to its stored power parameters
- A preset button will be dark gray with a white number when active , as shown for preset #1. (See figure 6.)
- Touch the S button followed by a numbered preset button to store the power parameters that are currently displayed

Adjustment Controls

- Adjust the highlighted orange setting (pulse energy, frequency or visual aiming beam intensity)
- Up arrow increases setting
- Down arrow decreases setting

Information Bar

- Displays simple instructions, system status and error messages

Pulse Width Toggle Button

- Toggles between Short and Long pulse width modes

Standby Button

- Yellow bar and icon signals system is in Standby mode
- Touch to return to Standby mode

Visual Aiming Beam Button


- Touch to toggle on or off
- The icon will be white when on and gray when off
- Adjust intensity with the adjustment controls while highlighted orange

Ready Button

- Green bar and icon signals system is in Ready mode
- Touch to enter Ready mode

Options Screen

The H-30 has many options that are displayed and set via the options screen. There are six option subscreens. The general layout of all the subscreens are identical, but the active screen is indicated by text in the information bar as well as the white text ("GEN" below) in the lower buttons. Touch the desired button to display the corresponding subscreen.

Options Screen Items (See figure 7)

Options Screen Tab

- Displays subscreens of options and the controls for adjustment
- Touch to enter options screen

Main Text Box

- Located in the upper left and displays relevant data

Adjustment Controls

- Control scrolling of main text box or active line (displayed in yellow text) within main text box
- Up arrow scrolls up
- Down arrow scrolls down

Information Bar

- Displays subscreen title and/or simple instructions

Enter Button

- Used to set various options as directed by instructions in the information bar
- Button text and icon will be white when active and gray when inactive

Laser Parameters

The following tables display available laser parameters; possible combination cells are marked with the average power (Power = Energy x Frequency) in watts for all available settings.

With a 150 µm fiber, the maximum power is 6 W

With a 200 µm fiber, the maximum power is 10 W

With a 273 µm fiber, the maximum power is 15 W

With a 365 µm fiber, the maximum power is 20 W

Combinations in Short pulse width mode:

Watts	5 Hz	8 Hz	10 Hz	12 Hz	15 Hz	20 Hz
0.5 J	2.5	4.0	5.0	6.0	7.5	10.0
0.6 J	3.0	4.8	6.0	7.2	9.0	
0.7 J	3.5	5.6	7.0	8.4	10.5	
0.8 J	4.0	6.4	8.0	9.6	12.0	
1.0 J	5.0	8.0	10.0	12.0	15.0	
1.2 J	6.0	9.6	12.0	14.4	18.0	
1.5 J	7.5	12.0	15.0	18.0	22.5	
2.0 J	10.0	16.0	20.0	24.0		
2.5 J	12.5	20.0	25.0	30.0		
3.0 J	15.0	24.0	30.0			
3.5 J	17.5					

Combinations in Long pulse width mode:

Watts	5 Hz	8 Hz	10 Hz	12 Hz	15 Hz	20 Hz
0.5 J	2.5	4.0	5.0	6.0	7.5	10.0
0.6 J	3.0	4.8	6.0	7.2	9.0	
0.7 J	3.5	5.6	7.0	8.4	10.5	
0.8 J	4.0	6.4	8.0	9.6	12.0	
1.0 J	5.0	8.0	10.0	12.0	15.0	
1.2 J	6.0	9.6	12.0	14.4	18.0	
1.5 J	7.5	12.0	15.0	18.0	22.5	
2.0 J	10.0	16.0	20.0	24.0		
2.5 J	12.5	20.0	25.0			
3.0 J	15.0					

Turning the System Off

To turn off the H-30, place the system into Standby mode, turn the key switch to the off (O) position and remove the key.

WARNING! If the system is turned off without entering Standby mode, the parameters will not be saved.

Protection Against Unauthorized Use

When the H-30 is not in use, protect it against unauthorized use by removing the key from the key switch.

Error Messages

There are several error messages that can be displayed either in the information bar on the main screen or in the error log on the options screen.

The following are possible types of messages:

- **Simmer failure:** The flash lamp doesn't simmer. Call for service.
- **Shutter failure:** The shutter isn't functioning properly. Call for service.
- **Over temperature:** The temperature of the cooling water is too high. Call for service.
- **Low coolant flow:** If water level is low, fill the reservoir with distilled water.

- **No fiber attached:** The fiber is not connected or is connected incorrectly. Unplug and then plug in the fiber correctly.
- **Remote open:** The remote connector, near the foot switch, is not connected or is connected incorrectly. This error warning is displayed when the system is in Ready mode. Connect remote.
- **Pedal not connected:** The foot switch connector is not connected or is connected incorrectly. This error warning is displayed when the system is in Ready mode. Connect foot switch.
- **Energy < 80%:** This warning appears when the output energy is less than 80% of the set value. Call for service.
- **Energy > 120%:** This warning appears when the output energy is greater than 120% of the set value. Call for service.

Chapter 6

Maintenance

General Maintenance

With careful use in normal operational conditions, the supplier recommends an overall system check every 12 months. Intense use, a dusty or dirty operational environment, or frequent movement requires more frequent maintenance. Verify the alignment and calibration using the instructions in Chapter 2.

WARNING! Only trained and qualified personnel should provide service for the H-30. Service by unauthorized personnel will void the warranty.

Cooling Water Level

Refill the water tank if the level of the water is too low.

Check the water level at least every six months.

Power Cord

The H-30 has a 15 ft. (5 m) power cord mechanically affixed to it.

The power cord may deteriorate over time because of general use and environmental conditions. Periodically check the status of the power cord and contact Customer Relations if a problem exists.

Labels

It is up to the user to maintain the safety labels and ensure that they are in good condition.

Replacing the Blast Shield

The blast shield is a replaceable optical window placed between the output focusing lens and the fiber bulkhead connector. Its purpose is to protect the output focusing lens from any debris in the case of fiber damage.

The blast shield is located on a bracket that is inserted from the top of the output focusing lens mount, just behind the laser aperture/fiber bulkhead connector.

Required Equipment: #2 Phillips screwdriver

Estimated Replacement Time: 2 minutes

Procedure

1. Unscrew and remove the decorative strip at the top of the rear panel and set aside. (See screw locations above.)



2. Completely loosen the locking screw on the bayonet latch. The locking screw is a captive screw, so it will not fall out when completely loosened.

3. Pull the bayonet latch back to open/disengage.

4. Lift the top panel (tabletop) and set the holding bar to hold the top panel up.

5. Unscrew the blast shield and slide it up, out of its mounting slot. The thumbscrew attached to the blast shield is a captive screw, so it will not fall off of the blast shield.
6. The blast shield can be inspected (for soiling or debris), and if damaged, should be replaced.
7. After replacing the blast shield, close the top panel. The blast shield can only be inserted and the thumbscrew attached in one orientation.
8. Engage the bayonet latch by pushing it in.
9. Loosely tighten the bayonet latch locking screw.
10. Reattach the decorative strip at the top of the rear panel.

Troubleshooting

Trouble	Possible Cause	Solution
Key switch does nothing	<ul style="list-style-type: none"> Power cord not connected Main line breaker off Emergency stop button pressed 	<ul style="list-style-type: none"> Connect power cord Check MAINS breaker Release the emergency stop button
No correspondence between displayed energy and effects on skin	<ul style="list-style-type: none"> Damaged fiber Damage blast shield 	<ul style="list-style-type: none"> Change fiber Change blast shield
Simmer failure	<ul style="list-style-type: none"> Flash lamp not functional 	<ul style="list-style-type: none"> Call for service

Error Messages

Error Message	Possible Cause	Solution
Low coolant flow	<ul style="list-style-type: none"> Not enough water 	<ul style="list-style-type: none"> Fill water reservoir
Over temperature	<ul style="list-style-type: none"> Room temperature too high Airflow is obstructed Cooling system component failure 	<ul style="list-style-type: none"> Lower the room temperature Ensure nothing is obstructing the fan Call for service
Shutter failure	<ul style="list-style-type: none"> Component failure 	<ul style="list-style-type: none"> Call for service
No fiber attached	<ul style="list-style-type: none"> Incorrect installation of the fiber Wrong fiber installed 	<ul style="list-style-type: none"> Check the proper connection of the fiber Install suitable H-30 fiber
Remote open	<ul style="list-style-type: none"> The remote interlock is not connected or is opened The remote interlock is damaged 	<ul style="list-style-type: none"> Check the remote interlock connection Call for service
Pedal not connected	<ul style="list-style-type: none"> The foot switch is not connected The foot switch is damaged 	<ul style="list-style-type: none"> Connect the foot switch properly Call for service
No energy	<ul style="list-style-type: none"> Optical component damage 	<ul style="list-style-type: none"> Call for service

Error Message	Possible Cause	Solution
Energy < 80%*	<ul style="list-style-type: none"> Optical component damage 	<ul style="list-style-type: none"> If the energy error persists, call for service*
Energy > 120%*	<ul style="list-style-type: none"> Incorrect energy calibration 	<ul style="list-style-type: none"> If the energy error persists, call for service*
Energy > 200%	<ul style="list-style-type: none"> Incorrect energy calibration 	<ul style="list-style-type: none"> Call for service
Maintenance due	<ul style="list-style-type: none"> Preventative maintenance is due 	<ul style="list-style-type: none"> Call for service*
MPS charge error	<ul style="list-style-type: none"> Main power supply malfunction 	<ul style="list-style-type: none"> Call for service
MPS over temp	<ul style="list-style-type: none"> Main power supply malfunction 	<ul style="list-style-type: none"> Call for service
MPS overload	<ul style="list-style-type: none"> Main power supply malfunction 	<ul style="list-style-type: none"> Call for service

*Minor error. Schedule service but ok to continue using laser.

Chapter 7

Warranty Information

[Return to Table of Contents](#)

Warranty Statement

Cook Urological Incorporated (“Company”) warrants to the purchaser of this H-30 Holmium Laser (the “Product”) that, at the time of manufacture, the Product was manufactured and tested in accordance with Good Manufacturing Practices specified by the United States Food and Drug Administration and was true to label. In the unlikely event of a Product failure due to defects in material or workmanship within the warranty period, the Product will be repaired or replaced, at Cook’s option, free of charge.

Because the Company has no control over the conditions under which this Product is used, diagnosis of the patient, anatomical or biological differences among patients, method of use or administration, or handling of the Product after it leaves the Company’s possession, the Company does not warrant a good effect or against an ill effect following the Product’s use. The foregoing warranty is exclusive and in lieu of all other warranties whether written, oral or implied. There are no warranties of merchantability or fitness for a particular purpose. No representative of the Company may change any of the foregoing, and the purchaser accepts the Product subject to all terms hereof. The Company reserves the right to change or discontinue this Product without notice.

Limits and Exclusions

The warranty does not cover consumables such as fiberoptic delivery systems, blast shields and output lenses, or any claims related to lost profits, loss of use, or any other indirect or consequential damages of any kind.

The warranty does not apply in the event of misuse, negligence or accidental damage of Product. Additionally, the following events are specifically excluded from warranty coverage: system damage due to plugging the Product into an improper electrical outlet, enclosure damage due to mechanical shock or exposure to high temperatures, laser cavity damage due to use of inappropriate or insufficient coolant and/or damage from exposure to freezing conditions.

The warranty will be null and void under any of the following conditions: an unauthorized service person has attempted to service the Product, unauthorized modifications have been made to the Product, use of uncertified fiberoptic accessories with Product, use of uncertified service parts with Product, and/or the Product’s optical system, pulse counter or any electronic safety circuit have been tampered with.

Warranty services will be provided only to the original purchaser at the location where the Product was originally installed. The warranty may not be transferred without the prior written consent of the Company.

Warranty Period

This limited warranty begins on the date the Product ships from the Company's facilities and ends on the same day 12 months later. For replacement Products provided during a warranty period, the remaining portion of the original warranty period shall apply to the replacement Product.

Warranty Procedure

Many Product issues can be resolved remotely. If you purchased your Product from a distributor, please contact your distributor for technical assistance, repair and replacement assistance.

If you purchased the Product from Cook, before returning a Product to Cook for service or repair, please contact Cook to speak with a service representative. If an issue cannot be resolved remotely, Cook may, at its discretion, repair or replace the defective Product or component at customer's location or at a Cook-designated repair facility. If the repairs will be made at a Cook-designated repair facility, purchaser is responsible for the cost of shipping of the Product to Cook's designated repair facilities, and Cook will pay return shipping back to the purchaser. Customer is responsible for appropriate packaging, insurance and risk of damage during transit to Cook's designated repair facility.

Warranty does not guarantee repair response times.

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VISIBLE AND INVISIBLE LASER RADIATION

Avoid eye or skin exposure to direct or scattered radiation.
 Class 4 laser product.

Maximum Power	30W
Pulse duration	150-800 μ s
Wavelengths	2100 nm
Aiming beam (Class 3R)	< 5 mW @ 532nm

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Standard IEC EN 60825-1: 2007

LASER RADIATION
 AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION

CLASS 4 LASER PRODUCT

CAUTION – CLASS 4 LASER RADIATION
 WHEN OPEN AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION

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ATTENTION – HIGH VOLTAGE
 THE CAPACITORS CAN BE CHARGED AT A POTENTIALLY DANGEROUS ELECTRIC VOLTAGE BEFORE REMOVING THE PROTECTIONS:

- SWITCH OFF POWER SUPPLY
- WAIT FOR 5 MINUTES
- CHECK FOR RESIDUAL VOLTAGE ON CAPACITORS

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Bloomington, IN 47404, U.S.A
COOK IRELAND LTD.