

## PERICLES, Proximal Femoral Nail

### Product Information

#### Nails

REF. NO.	Product Description	Specification
070890170	PERICLES Proximal Femoral Nail	φ9.0x170mm
070890200	PERICLES Proximal Femoral Nail	φ9.0x200mm
070890240	PERICLES Proximal Femoral Nail	φ9.0x240mm
070900170	PERICLES Proximal Femoral Nail	φ10.0x170mm
070900200	PERICLES Proximal Femoral Nail	φ10.0x200mm
070900240	PERICLES Proximal Femoral Nail	φ10.0x240mm
070910170	PERICLES Proximal Femoral Nail	φ11.0x170mm
070910200	PERICLES Proximal Femoral Nail	φ11.0x200mm
070910240	PERICLES Proximal Femoral Nail	φ11.0x240mm
070920170	PERICLES Proximal Femoral Nail	φ12.0x170mm
070920200	PERICLES Proximal Femoral Nail	φ12.0x200mm
070920240	PERICLES Proximal Femoral Nail	φ12.0x240mm
071111300	PERICLES Proximal Femoral Nail, Long, Left	φ9.0x300mm
071111340	PERICLES Proximal Femoral Nail, Long, Left	φ9.0x340mm
071111380	PERICLES Proximal Femoral Nail, Long, Left	φ9.0x380mm
071111420	PERICLES Proximal Femoral Nail, Long, Left	φ9.0x420mm
071112300	PERICLES Proximal Femoral Nail, Long, Right	φ9.0x300mm
071112340	PERICLES Proximal Femoral Nail, Long, Right	φ9.0x340mm
071112380	PERICLES Proximal Femoral Nail, Long, Right	φ9.0x380mm
071112420	PERICLES Proximal Femoral Nail, Long, Right	φ9.0x420mm
071121300	PERICLES Proximal Femoral Nail, Long, Left	φ10.0x300mm
071121340	PERICLES Proximal Femoral Nail, Long, Left	φ10.0x340mm
071121380	PERICLES Proximal Femoral Nail, Long, Left	φ10.0x380mm
071121420	PERICLES Proximal Femoral Nail, Long, Left	φ10.0x420mm
071122300	PERICLES Proximal Femoral Nail, Long, Right	φ10.0x300mm
071122340	PERICLES Proximal Femoral Nail, Long, Right	φ10.0x340mm
071122380	PERICLES Proximal Femoral Nail, Long, Right	φ10.0x380mm
071122420	PERICLES Proximal Femoral Nail, Long, Right	φ10.0x420mm
071131300	PERICLES Proximal Femoral Nail, Long, Left	φ11.0x300mm
071131340	PERICLES Proximal Femoral Nail, Long, Left	φ11.0x340mm
071131380	PERICLES Proximal Femoral Nail, Long, Left	φ11.0x380mm
071131420	PERICLES Proximal Femoral Nail, Long, Left	φ11.0x420mm
071132300	PERICLES Proximal Femoral Nail, Long, Right	φ11.0x300mm
071132340	PERICLES Proximal Femoral Nail, Long, Right	φ11.0x340mm
071132380	PERICLES Proximal Femoral Nail, Long, Right	φ11.0x380mm
071132420	PERICLES Proximal Femoral Nail, Long, Right	φ11.0x420mm
071141300	PERICLES Proximal Femoral Nail, Long, Left	φ12.0x300mm
071141340	PERICLES Proximal Femoral Nail, Long, Left	φ12.0x340mm
071141380	PERICLES Proximal Femoral Nail, Long, Left	φ12.0x380mm
071141420	PERICLES Proximal Femoral Nail, Long, Left	φ12.0x420mm
071142300	PERICLES Proximal Femoral Nail, Long, Right	φ12.0x300mm
071142340	PERICLES Proximal Femoral Nail, Long, Right	φ12.0x340mm
071142380	PERICLES Proximal Femoral Nail, Long, Right	φ12.0x380mm
071142420	PERICLES Proximal Femoral Nail, Long, Right	φ12.0x420mm

#### Others

REF. NO.	Product Description	Length(mm)
070370075-070370120	PERICLES Spiral Blade	75-120 (5mm increments)
070120025-070120085	4.9mm Locking Screw	25-85 (5mm increments)
070420000	PERICLES End Cap	
070421000	PERICLES End Cap, Long	





## PERICLES, Proximal Femoral Nail



7s Medical AG  
Münsterstrasse 1i  
CH- 6214 Schenkon  
Switzerland

Phone +41 41 925 66 60  
info@7smedical.com  
www.7smedical.com

**PERICLES, Proximal Femoral Nail**






Increased blade-bone interface provides better stability in anchoring and antirotation.

Lemon shaft design provides antirotation mechanism at the blade – nail junction.

Distal flutes facilitates the insertion of the nail reducing the stress on the bone at the tip of the nail.

Allows static and dynamic distal locking options for clinical choice.



PERICLES diameter 17mm.

Lateral flat cut design for easier insertion and lower pressure on the lateral cortex.

Lateral Medial 6 degree angulation allows an entrance point at the top of the great trochanter.

**Indications & Contraindications**

PERICLES Proximal Femur Nail

**Short (170mm – 240mm)**

Indications

- Pertrochanteric fractures
- Intertrochanteric fractures
- High subtrochanteric fractures

Contraindications

- Low subtrochanteric fractures
- Femoral shaft fractures
- Isolated or combined medial femoral neck fractures

**Long (300mm – 420mm)**

Indications

- Low and extended subtrochanteric fractures
- Ipsilateral trochanteric fractures
- Combination fractures (in the proximal femur)
- Pathological fractures

Contraindications

- Isolated femoral shaft fractures
- Isolated or combined medial femoral neck fractures



4.9mm Locking Screw

**Warning**

This flyer is just for understanding the specific product features. For clinical usage, please refer to the surgical guide. Instruction by experienced surgeon is highly recommended.