


COOK[®]

MEDICAL


TURP/TURBT

RESECTION ELECTRODES




RESECTOSCOPE ELECTRODES

Cook offers a full line of resectoscope electrodes for urological procedures including coagulating, cutting, vaporization, roller ball, roller barrel, and knife. Our electrodes are compatible with most major resectoscopes.





Cutting Loop

Made with the same quality synonymous with Cook. Our tungsten wire cutting loops are hand-formed.




Coagulating Perforated Resector

Produces the chip size of a cutting loop but provides a larger zone of coagulation which reduces intraoperative bleeding.




Roller Ball and Roller Barrel

Used for coagulation. Cook offers 3 & 5mm roller balls and 3mm roller barrels.



Collings Knife

For access to all areas. Commonly used for ureteral meatotomy and transurethral incision of the prostate (TUIP).



Dimpled Vaporization Roller

Provides a smoother area of vaporization than a grooved roller while simultaneously ablating and coagulating tissue.

All products supplied 5 per box.

	Power Type	Valleylab Ref
Cutting Loop	Cut, Coag, Blend	TURP
Perforated Resector	Cut	TURVP
Roller Ball/Barrel	Coag	TURVP
Dimpled Roller	Cut	TURVP
Knife	Cut, Coag, Blend	TURP

Cook Description	Cook Order #	Microvasive Description	Microvasive Order #	Bard Description	Bard Order #	ACMI Description	ACMI Order #
Cutting Loop 24 Fr, 0.012 Wire	CL-2412-AC	Cutting Loop (24 Fr)	M0068802131	C-Max Cutting Loop Electrode (24 Fr)	355301	Right Angle Cutting Loop	MLE-24-012
Cutting Loop 26 Fr, 0.012 Wire	CL-2612-AC	Cutting Loop (26 Fr)	M0068802141	C-Max Cutting Loop Electrode (26 Fr)	355311	Right Angle Cutting Loop	MLE-26-012
Cutting Loop 28 Fr, 0.012 Wire	CL-2812-AC	Cutting Loop (28 Fr)	M0068802191	C-Max Cutting Loop Electrode (28 Fr)	355321	Right Angle Cutting Loop	MLE-28-012
Cutting Loop 24 Fr, 0.015 Wire	CL-2415-AC					Right Angle Cutting Loop	MLE-24-015
Cutting Loop 26 Fr, 0.015 Wire	CL-2615-AC					Right Angle Cutting Loop	MLE-26-015
Cutting Loop 28 Fr, 0.015 Wire	CL-2815-AC					Right Angle Cutting Loop	MLE-28-015
Coagulating Perforated Resector 24 Fr	PR-2400-AC	Titanium Wedge (24 Fr)	M0068803111	Max-Blade Coagulating Resector (24 Fr)	355201	Coagulating Electrode	CE-24
						Grooved VaporTome Vaporization Electrode	VE-LG
Coagulating Perforated Resector 26 Fr	PR-2600-AC	Titanium Wedge (26 Fr)	M0068803121	Max-Blade Coagulating Resector (26 Fr)	355211	Coagulating Electrode	CE-26
						Grooved VaporTome Vaporization Electrode	VE-LG
Coagulating Perforated Resector 28 Fr	PR-2800-AC			Max-Blade Coagulating Resector (28 Fr)	355221	Coagulating Electrode	CE-28
						Grooved VaporTome Vaporization Electrode	VE-LG
3 mm Roller Ball 24 Fr (Fits 24, 26, 28 Fr)	RB-243M-AC	3 mm Rollerball Electrode (24 Fr)	M0068802151			Roller Ball Electrode	RE RE-28
		3 mm Rollerball Electrode (26 Fr)	M0068802161				
5 mm Roller Ball 24 Fr (Fits 24, 26, 28 Fr)	RB-245M-AC	5 mm Rollerball Electrode (24 Fr)	M0068804111				
		5 mm Rollerball Electrode (26 Fr)	M0068804121				
3 mm Roller Barrel 24 Fr (Fits 24, 26, 28 Fr)	RL-243M-AC					Roller Bar Electrode	RB
3mm Dimpled Vaporization Roller 24 Fr (Fits 24, 26, 28 Fr)	VR-243M-AC	Powerball Electrode (24 Fr)	M0068802111	V-Max Roller Bar Electrode (24 Fr)	355401	Grooved VaporTrome Vaporization Electrode	VE-B
						Fluted VaporTrome Vaporization Electrode	VE-F
						Spiked VaporTrome Vaporization Electrode	VE-S
Collings Knife 24 Fr (Fits 24, 26, 28 Fr)	CK-2400-AC	TUIP Knife (24 Fr)	M0068802171			PERKASH Right Angle Knife Electrode	KE-24 KE-26 KE-28
		Incision Blade (24 Fr)	M0068805111			PERKASH Right Angle Knife Electrode	RAK-24 RAK-26 RAK-28

	Power Type
	Cut, Coag, Blend
Valleylab Ref	Cut
TURP	Coag
TURVP	Cut
TURVP	Cut, Coag, Blend
Cutting Loop	
Perforated Resector	
Roller Ball/Barrel	
Dimpled Roller	
Knife	

Cook Description	Cook Order #	Microvasive Description	Microvasive Order #	Bard Description	Bard Order #	Olympus Description	Olympus Order #
Cutting Loop 24 Fr, 0.012 Wire, 30°	CL-2412-30-OL	Cutting Loop (24 Fr)	M0068802331			30° Cutting Loop	A22205C
Cutting Loop 26 Fr, 0.012 Wire, 30°	CL-2612-30-OL	Cutting Loop (26 Fr)	M0068802341			30° Cutting Loop	A22206C
Cutting Loop 28 Fr, 0.012 Wire, 30°	CL-2812-30-OL	Cutting Loop (28 Fr)	M0068802391			30° Cutting Loop	A22207C
Cutting Loop 24 Fr, 0.012 Wire, 12°	CL-2412-12-OL			C-Max Cutting Loop Electrode (24 Fr)	355303	12° Cutting Loop	A22201C
Cutting Loop 26 Fr, 0.012 Wire, 12°	CL-2612-12-OL			C-Max Cutting Loop Electrode (26 Fr)	355313	12° Cutting Loop	A22202C
Cutting Loop 28 Fr, 0.012 Wire, 12°	CL-2812-12-OL			C-Max Cutting Loop Electrode (28 Fr)	355323	12° Cutting Loop	A22203C
Coagulating Perforated Resector 24 Fr	PR-2400-OL	Titanium Wedge (24 Fr)	M0068803311	Max-Blade Coagulating Resector (24 Fr)	355203	Band Loop	A22223C A22221C
Coagulating Perforated Resector 26 Fr	PR-2600-OL	Titanium Wedge (26 Fr)	M0068803321	Max-Blade Coagulating Resector (26 Fr)	355213	Band Loop	A22224C A22222C
Coagulating Perforated Resector 28 Fr	PR-2800-OL			Max-Blade Coagulating Resector (28 Fr)	355223		
3 mm Roller Ball 24 Fr (Fits 24, 26, 28 Fr)	RB-243M-OL	3 mm Rollerball Electrode (24 Fr)	M0068802351			3 mm Roller Ball Electrode	A22251C
		3 mm Rollerball Electrode (26 Fr)	M0068802361				
5 mm Roller Ball 24 Fr (Fits 24, 26, 28 Fr)	RB-245M-OL	5 mm Rollerball Electrode (24 Fr)	M0068804311			Large Roller Ball Electrode (5 mm)	A22258C
		5 mm Rollerball Electrode (26 Fr)	M0068804321				
3 mm Roller Barrel 24 Fr (Fits 24, 26, 28 Fr)	RL-243M-OL					Cylinder Electrode (3 mm)	A 2196 A 4753
3 mm Dimpled Vaporization Roller 24 Fr (Fits 24, 26, 28 Fr)	VR-243M-OL	Powerball Electrode (24 Fr)	M0068802311	V-Max Roller Bar Electrode (24 Fr)	355403	Grooved Cylinder Electrode	A22257C A22257D
Collings Knife 24 Fr (Fits 24, 26, 28 Fr)	CK-2400-OL	TUIP Knife (24 Fr)	M0068802371			Knife Electrode	A22265C
		Incision Blade (24 Fr)	M0068805311				

	Power Type	Valleylab Ref	
Cutting Loop	Cut, Coag, Blend	TURP	
Perforated Resector	Cut	TURVP	
Roller Ball/Barrel	Coag	TURVP	
Dimpled Roller	Cut	TURVP	
Knife	Cut, Coag, Blend	TURP	

Cook Description	Cook Order #	Microvase Description	Microvase Order #	Bard Description	Bard Order #	Storz Description	Storz Order #
Cutting Loop 24 Fr, 0.012 Wire	CL-2412-KS1	Cutting Loop (24 Fr)	M0068802031	C-Max Cutting Loop Electrode (24 Fr)	355305	Angled Cutting Loop	27050 G
Cutting Loop 27 Fr, 0.012 Wire	CL-2712-KS1	Cutting Loop (27 Fr)	M0068802041	C-Max Cutting Loop Electrode (27 Fr)	355315	Angled Cutting Loop	27050 F
Cutting Loop 24 Fr, 0.014 Wire	CL-2414-KS1			C-Max Cutting Loop Electrode (24 Fr)	355325	Angled Cutting Loop	27050 G
Cutting Loop 27 Fr, 0.014 Wire	CL-2714-KS1			C-Max Cutting Loop Electrode (27 Fr)	355335	Angled Cutting Loop	27050 F
Coagulating Perforated Resector 24 Fr	PR-2400-KS1	Titanium Wedge (24 Fr)	M0068803011	Max-Blade Coagulating Resector (24 Fr)	355205	VaporCut Electrode	27050 WG 27050 SG 27050 GB
Coagulating Perforated Resector 27 Fr	PR-2700-KS1	Titanium Wedge (27 Fr)	M0068803021	Max-Blade Coagulating Resector (27 Fr)	355215	VaporCut Electrode	27050 WF 27050 SF 27050 GF
3 mm Roller Ball 24 Fr (Fits 24, 27 Fr)	RB-243M-KS1	3 mm Rollerball Electrode (24 Fr)	M0068802051			3 mm Coagulating Ball End Electrode	27050 N
		3 mm Rollerball Electrode (27 Fr)	M0068802061			3 mm Coagulating Ball End Electrode	27050 M
5 mm Roller Ball 24 Fr (Fits 24, 27 Fr)	RB-245M-KS1	5 mm Rollerball Electrode (24 Fr)	M0068804011			5 mm Coagulating Ball End Electrode	27050 NK
		5 mm Rollerball Electrode (27 Fr)	M0068804021			5 mm Coagulating Ball End Electrode	27050 MK
3 mm Roller Barrel 24 Fr (Fits 24, 27 Fr)	RL-243M-KS1					3 mm Coagulating Barrel End Electrode	27050 NX 27050 MX
						5 mm Coagulating Barrel End Electrode	27050 NW 27050 MW
3 mm Dimpled Vaporization Roller 24 Fr (Fits 24, 27 Fr)	VR-243M-KS1	Powerball Electrode (24 Fr)	M0068802011	V-Max Roller Bar Electrode (24 Fr)	355404		
Collings Knife 24 Fr (Fits 24, 27 Fr)	CK-2400-KS1	TUIP Knife (24 Fr)	M0068802071			Coagulating Electrode Pointed	27050 L
		Incision Blade (24 Fr)	M0068805011			Coagulating Electrode Pointed	27050 K

	Valley/lab Ref	Power Type
Cutting Loop	TURP	Cut, Coag, Blend
Perforated Resector	TURVP	Cut
Roller Ball/Barrel	TURVP	Coag
Dimpled Roller	TURVP	Cut
Knife	TURP	Cut, Coag, Blend

Cook Description	Cook Order #	Microvase Description	Microvase Order #	Bard Description	Bard Order #	Storz Description	Storz Order #
Cutting Loop 24 Fr, 0.014 Wire	CL-2414-KS2	Cutting Loop (24 Fr)	M0068802231	C-Max Cutting Loop Electrode (24 Fr)	355324	Angled Cutting Loop	27040 G
Cutting Loop 27 Fr, 0.014 Wire	CL-2714-KS2	Cutting Loop (27 Fr)	M0068802241	C-Max Cutting Loop Electrode (27 Fr)	355334	Angled Cutting Loop	27040 F
Coagulating Perforated Resector 24 Fr	PR-2400-KS2	Titanium Wedge (24 Fr)	M0068803211	Max-Blade Coagulating Resector (24 Fr)	355204	VaporCut Electrode	27040 WG 27040 SG 27040 GB
Coagulating Perforated Resector 27 Fr	PR-2700-KS2	Titanium Wedge (27 Fr)	M0068803221	Max-Blade Coagulating Resector (27 Fr)	355214	VaporCut Electrode	27040 WF 27040 SF 27040 GF
3 mm Roller Ball 24 Fr (Fits 24, 27 Fr)	RB-243M-KS2	3 mm Rollerball Electrode (24 Fr)	M0068802251			3 mm Coagulating Ball End Electrode	27040 N 27040 M
		3 mm Rollerball Electrode (27 Fr)	M0068802261				
5 mm Roller Ball 24 Fr (Fits 24, 27 Fr)	RB-245M-KS2	5 mm Rollerball Electrode (24 Fr)	M0068804211			5 mm Coagulating Ball End Electrode	27040 NK 27040 MK
		5 mm Rollerball Electrode (27 Fr)	M0068804221				
3 mm Roller Barrel 24 Fr (Fits 24, 27 Fr)	RL-243M-KS2					3 mm Coagulating Barrel End Electrode	27040 NX 27040 MX
						5 mm Coagulating Barrel End Electrode	27040 NW 27040 MW
3 mm Dimpled Vaporization Roller 24 Fr (Fits 24, 27 Fr)	VR-243M-KS2	Powerball Electrode (24 Fr)	M0068802211	V-Max Roller Bar Electrode (24 Fr)	355404		