

Contents – Catheters

Ureteral catheters

Dual Lumen Ureteral Access Catheter	2.2
Cone Tip Ureteral Catheter	2.3
Flexi-Tip® Ureteral Catheter – Closed End.	2.4
Open-End Flexi-Tip® Ureteral Catheter	2.5
Angled Tip Ureteral Catheter.	2.6
Open-End Ureteral Catheter	2.7
Open-End Ureteral Catheter – Sof-Flex®	2.8
Whistle Tip Ureteral Catheter.	2.9

Pigtail ureteral catheters

Pigtail Ureteral Catheter Set.	2.10
--	------

Renal access catheters

Renal Access Cobra Catheter	2.11
---------------------------------------	------



Dual Lumen Ureteral Access Catheter


Used for the injection of contrast medium and/or for safety wire placement during percutaneous and transurethral access. The catheter is available with an AQ® coating. This is a microthin, hydrophilic polymer coating that, when activated, attracts and holds water and other liquids to the catheter to create a low-friction surface.



Order Number	Reference Part Number	Fr	Length cm	Placement Lumen Diameter inch	Injection Lumen Diameter inch	Tip Configuration
Standard						
G17013	022510	5.0/10.0	50	0.050	0.050	
G18845	022610-24	6.0/10.0	24	0.050	0.050	Flexi-Tip®
G17145	022610	6.0/10.0	50	0.050	0.050	Flexi-Tip®
With AQ Hydrophilic Coating						
G17560	AQ-022510	5.0/10.0	50	0.050	0.050	
G17323	AQ-022610	6.0/10.0	50	0.050	0.050	Flexi-Tip®
With Radiopaque Marker						
G17330	022610-RM	6.0/10.0	50	0.050	0.050	Flexi-Tip®


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Cone Tip Ureteral Catheter

Used for retrograde pyelogram. The tip occludes the ureteral orifice.



Order Number	Reference Part Number	Fr	Length cm	Tip Fr	Accepts Wire Guide Diameter inch
With Closed End					
G14663	024605	4.8	70	8.0	-
G14665	024606	6.0	70	10.0	-
With Open End					
G14717	024605-S1	4.8	70	8.0	0.035

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Flexi-Tip® Ureteral Catheter

Closed End

Indicated for access and catheterisation of the urinary tract, including the following applications:


- Delivery of contrast media
- Drainage of fluids from the urinary tract
- Delivery of irrigation fluids to the urinary tract
- Navigation of a tortuous ureter



Order Number	Reference Part Number	Fr	Length cm	Quantity
Right				
G14014	021103	3.0	70	10
G14015	021104	4.0	70	10
G14016	021105	5.0	70	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Open-End Flexi-Tip® Ureteral Catheter

Indicated for access and catheterisation of the urinary tract, including the following applications:

- Delivery of contrast media
- Drainage of fluids from the urinary tract
- Delivery of irrigation fluids to the urinary tract
- Navigation of a tortuous ureter
- Access, advancement, or exchange of wire guides



Order Number	Reference Part Number	Fr	Length cm	Accepts Instrumentation Diameter inch/mm	Comment	Quantity
Standard						
G15040	021304	4.0	70	0.028/0.71		10
G14521	021305	5.0	70	0.038/0.97		10
G14808	021306	6.0	70	0.038/0.97		10
G14809	021307	7.0	70	0.038/0.97		10
Modified						
G15889	021305-S13	5.0	70	0.038/0.97	2 sideports, bevelled	1
G17051	021305-S18	5.0	10	0.038/0.97	10 cm in length	1

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.




Angled Tip Ureteral Catheter

Indicated for access and catheterisation of the urinary tract, including the following applications:

- Delivery of contrast media
- Navigation of a tortuous ureter
- Access, advancement, or exchange of wire guides


The catheter comes with a 0.035 inch diameter, 145 cm long stainless steel wire guide. The tip has a 20° angle.



Order Number	Reference Part Number	Fr	Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Tip Configuration
20° Angled Tip						
G14598	023105	5.0	70	0.035	145	tapered
G14587	023106	6.0	70	0.035	145	tapered

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.


Please see product risk information in the IFU at cookmedical.eu.



Open-End Ureteral Catheter

Indicated for access and catheterisation of the urinary tract, including the following applications:


- Delivery of contrast media
- Drainage of fluids from the urinary tract
- Delivery of irrigation fluids to the urinary tract
- Navigation of a tortuous ureter
- Access, advancement, or exchange of wire guides



Order Number	Reference Part Number	Fr	Length cm	Accepts Instrumentation inch/mm	Quantity
G14818	020013	3.0	70	0.018/0.46	1
G14429	020014	4.0	70	0.025/0.64	10
G14274	020015	5.0	70	0.038/0.97	10
G14430	020116	6.0	70	0.038/0.97	10
G14415	020017	7.0	70	0.038/0.97	10
G14456	020018	8.0	70	0.038/0.97	10
Modified					
G16193	020015-S13	5.0	120	0.038/0.97	1

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.


Please see product risk information in the IFU at cookmedical.eu.



Open-End Ureteral Catheter Sof-Flex®

Indicated for access and catheterisation of the urinary tract, including the following applications:


- Delivery of contrast media
- Drainage of fluids from the urinary tract
- Delivery of irrigation fluids to the urinary tract
- Navigation of a tortuous ureter
- Access, advancement, or exchange of wire guides



Order Number	Reference Part Number	Fr	Length cm	Accepts Instrumentation inch/mm	Comment
G15946	020025-C4	4.0	70	0.025/0.64	
G15943	020038-C5	5.0	70	0.038/0.97	
G15944	020040-C5	5.0	70	0.040/1.02	packaged coiled
G15942	020040-C6	6.0	70	0.040/1.02	packaged coiled

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.


Please see product risk information in the IFU at cookmedical.eu.



Whistle Tip Ureteral Catheter

Indicated for access and catheterisation of the urinary tract, including the following applications:


- Delivery of contrast media
- Drainage of fluids from the urinary tract
- Delivery of irrigation fluids to the urinary tract
- Navigation of a tortuous ureter



Order Number	Reference Part Number	Fr	Length cm	Quantity	Color
Left					
G14036	022203	3.0	70	10	Green
G14037	022204	4.0	70	10	Blue
G14038	022205	5.0	70	10	Yellow
Right					
G14028	022103	3.0	70	10	Green
G14029	022104	4.0	70	10	Blue
G14030	022105	5.0	70	10	Yellow
G14031	022106	6.0	70	10	Green


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Pigtail Ureteral Catheter Set

Indicated for drainage or irrigation in the urinary tract as well as access, advancement, or exchange of wire guides. The catheter has sideports in the pigtail and comes with a stainless steel wire guide.



Order Number	Reference Part Number	Fr	Length cm	Wire Guide Diameter inch	Wire Guide Length cm
Sideports in Pigtail Only					
Left					
G14148	025205	5.0	70	0.038	145
G14149	025206	6.0	70	0.038	145
G14150	025207	7.0	70	0.045	145
Right					
G14144	025105	5.0	70	0.038	145
G14145	025106	6.0	70	0.038	145
G14146	025107	7.0	70	0.045	145
G14803	025108	8.0	70	0.045	145


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Renal Access Cobra Catheter

Used percutaneously to direct a 0.038 inch (0.97 mm) diameter wire guide to a specific renal area for the delivery of contrast media and for the access, advancement, or exchange of wire guides. The catheter features Torcon® catheter material with an inner stainless steel braid.



Order Number	Reference Part Number	Fr	Length cm
G14408	023806	6.0	65

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

