

Contents – Dilation and access

Ureteral access sheaths

Flexor® Parallel™ Rapid Release™ Ureteral Access Sheath . .	4.2
Flexor® Ureteral Access Sheath	4.3
Safety Wire Guide Introducer Set	4.4

Access catheters

Desilets-Hoffman Introducer Set	4.5
Cook 810 Set™	4.6
Peel-Away® Introducer Set	4.7

Ureteral dilators


AQ® Long Taper Ureteral Dilator	4.8
Ureteral Dilator Set	4.9

Renal dilators

Amplatz Renal Dilator Set	4.10
Amplatz Renal Introducer	4.11
Amplatz Sheath – Clear with Radiopaque Stripe	4.12
Amplatz Sheath – Opaque	4.13
Fascial Dilator Set	4.14


Urethral dilators

Integral Filiiform Urethral Dilator	4.15
Tunner Paediatric Urethral Dilator Set	4.16
S~Curve™ Urethral Dilator Set	4.17
Urethral Dilator Set	4.18



Flexor® Parallel™ Rapid Release™ Ureteral Access Sheath

Used to establish a conduit during endoscopic urological procedures, facilitating the passage of endoscopes and other instruments into the urinary tract. The Flexor Parallel offers two options for placement: the traditional Flexor placement or the Rapid Release technique. The working channel protects the ureter during repeated instrument exchanges. The sheath also protects delicate instruments and smaller flexible scopes from damage. The sheath has a hydrophilic coating to create a low-friction surface as the ureter is dilated.




Order Number	Reference Part Number	Sheath ID Fr	Length cm
G56428	FUS-095020-P	9.5	20
G56429	FUS-095028-P	9.5	28
G56430	FUS-095035-P	9.5	35
G56431	FUS-095045-P	9.5	45
G56432	FUS-095055-P	9.5	55
G34871	FUS-107020-P	10.7	20
G34872	FUS-107028-P	10.7	28
G34873	FUS-107035-P	10.7	35
G34874	FUS-107045-P	10.7	45
G34875	FUS-107055-P	10.7	55

Order Number	Reference Part Number	Sheath ID Fr	Length cm
G56433	FUS-120020-P	12.0	20
G56434	FUS-120028-P	12.0	28
G56435	FUS-120035-P	12.0	35
G56436	FUS-120045-P	12.0	45
G56438	FUS-120055-P	12.0	55
G56439	FUS-140020-P	14.0	20
G56440	FUS-140028-P	14.0	28
G56441	FUS-140035-P	14.0	35
G56442	FUS-140045-P	14.0	45
G56443	FUS-140055-P	14.0	55


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Flexor® Ureteral Access Sheath

Used to establish a conduit during endoscopic urological procedures, facilitating the passage of endoscopes and other instruments into the urinary tract. The working channel protects the ureter during repeated instrument exchanges. The sheath also protects delicate instruments and smaller flexible scopes from damage. The sheath has a hydrophilic coating to create a low-friction surface as the ureter is dilated.



Order Number	Reference Part Number	Sheath ID Fr	Length cm
G32853	FUS-095013*	9.5	13
G19177	FUS-095020	9.5	20
G19176	FUS-095028	9.5	28
G19174	FUS-095035	9.5	35
G46167	FUS-095045	9.5	45
G19175	FUS-095055	9.5	55
G34876	FUS-107013*	10.7	13
G34877	FUS-107020	10.7	20
G34878	FUS-107028	10.7	28
G34879	FUS-107035	10.7	35
G34880	FUS-107045	10.7	45
G34881	FUS-107055	10.7	55

Order Number	Reference Part Number	Sheath ID Fr	Length cm
G32852	FUS-120013*	12.0	13
G19173	FUS-120020	12.0	20
G19172	FUS-120028	12.0	28
G19168	FUS-120035	12.0	35
G46168	FUS-120045	12.0	45
G19169	FUS-120055	12.0	55
G30391	FUS-140020	14.0	20
G30393	FUS-140028	14.0	28
G30389	FUS-140035	14.0	35
G46169	FUS-140045	14.0	45
G30390	FUS-140055	14.0	55


* For paediatric use in children who are two years old or older.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

Safety Wire Guide Introducer Set

Used for antegrade placement of a safety wire guide prior to endourological procedures.




Order Number	Reference Part Number	Sheath Material	Sheath ID/OD Fr	Sheath Length cm	Obturator Length cm
Sheath and Obturator					
G14306	074112	polythene	8/12	22.2	25.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.



Desilets-Hoffman Introducer Set

Used for retrograde or antegrade ureteral access and subsequent ureteral stent placement. The "RB" version has a radiopaque band on the delivery sheath to aid in fluoroscopic visualisation during placement and increase visibility.




Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm	Sheath ID Fr	Sheath OD Fr	Sheath Length cm
G18063	DHIS-080030-DIL70	8.0	70	8.0	10.0	30
With Radiopaque Band						
G18089	DHIS-080030-RB-DIL88	8.0	88	8.0	10.0	30

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

Cook 810 Set™

Used for retrograde or antegrade ureteral access and to facilitate wire guide exchange.



Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm	Sheath ID Fr	Sheath Length cm
G15010	030000	8.0	88	8.5	50
G15079	030000-S1	8.0	88	8.5	30


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Peel-Away® Introducer Set

Used for catheter or instrument introduction. The sheath peels away and can be removed after catheter placement. The maximum diameter of the catheter or instrument to be introduced should be measured to ensure its passage through the sheath.



Order Number	Reference Part Number	Sheath ID Fr	Dilator Length cm	Sheath Length cm
Set				
G14440	073909	9.0	37	32
G14441	073910	10.0	37	32
G14442	073912	12.0	37	32

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



AQ® Long Taper Ureteral Dilator

Used for dilation of the ureter prior to ureteroscopy and/or stone manipulation. The long taper eases introduction through strictured or tortuous areas. AQ is a microthin, hydrophilic polymer coating that, when activated, attracts and holds water and other liquids to the dilator. This creates a low-friction surface.




Order Number	Reference Part Number	Distal Fr	Proximal Fr	Length cm
G17636	AQ-076012-TP	6.0	12.0	60

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Ureteral Dilator Set


Used for dilation of the ureter prior to ureteroscopy and/or stone manipulation. AQ is a microthin, hydrophilic polymer coating that, when activated, attracts and holds water and other liquids to the dilator. This creates a low-friction surface.



Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm
Set			
Polythene with AQ Hydrophilic Coating			
G17310	AQ-076000	6.0–18.0	60
Polythene			
G14297	076000	6.0–18.0	60
PTFE			
G14574	076600	6.0–18.0	60


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Amplatz Renal Dilator Set

Used for nephrostomy tract dilation.



Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm	Sheath Length cm
Set				
G14292	075000	6.0–30.0	30	16
Modified				
G14485	075000-S1	6.0–30.0	30	20
With AQ Hydrophilic Coating				
G17322	AQ-075000	6.0–30.0	30	16
Guide Catheter Only				
G14305	075001	8.0	84	–


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Amplatz Renal Introducer


Used for nephrostomy tract dilation and the placement of a sheath.



Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm	Sheath Length cm
G18013	ARI-180016	18.0	30	16
G18001	ARI-240016	24.0	30	16
G18004	ARI-300016	30.0	30	16
G18016	ARI-300025	30.0	30	25

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Amplatz Sheath

Clear with Radiopaque Stripe


Used as a working channel to maintain a previously established nephrostomy tract.



Order Number	Reference Part Number	Sheath ID Fr	Length cm
G46030	CAS-180017	18.0	17
G46031	CAS-240017	24.0	17
G46032	CAS-300017	30.0	17

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Amplatz Sheath

Opaque

Used as a working channel to maintain a previously established nephrostomy tract.




Order Number	Reference Part Number	Sheath ID Fr	Length cm
G17984	VS-180016-ARI	18.0	16
G17964	VS-240016-ARI	24.0	16
G17967	VS-300016-ARI	30.0	16
G17976	VS-300030-ARI	30.0	30

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

Fascial Dilator Set


Used for fascial dilation.



Order Number	Reference Part Number	Dilator OD Fr	Length cm
Set			
PTFE			
G14572	073600	6.0–18.0	37
Polythene			
G14744	072000	6.0–24.0	20
G14214	073700	8.0–24.0	37
Dilator only			
G14133	072006	6.0	20
G14135	072008	8.0	20
G14137	072010	10.0	20
G14138	072012	12.0	20
G14139	072014	14.0	20
G14140	072016	16.0	20
G14141	072018	18.0	20


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Integral Filiform Urethral Dilator

Used for traversing and dilating urethral strictures. The filiform tip replaces the need for separate filiform and followers to dilate the urethra.



Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm	Filiform Tip Fr	Filiform Tip Length cm
G17909	FFD-180032	18.0	32	7.0	6


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

Tunner Paediatric Urethral Dilator Set




Used for paediatric urethral dilation. The radiopaque PTFE dilators have a closed, rounded, 20° angled tip.



Order Number	Reference Part Number	Dilator OD Fr	Length cm
G15262	073300	8.0–14.0	14


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



S~Curve™ Urethral Dilator Set


Used for dilation of male urethral strictures and vesical neck contractures. The S~Curve dilators support the natural curvature of the male urethra. The dilators have an AQ® hydrophilic coating that attracts and holds liquids to the dilator. This creates a low-friction, slippery surface to ease placement.



Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm	Comment
Set				
G32789	073701-CD	8.0–20.0	37	with sideport
Dilator Only				
G46249	073722-CD	22.0	37	with sideport
G46250	073724-CD	24.0	37	with sideport
G46251	073726-CD	26.0	37	with sideport


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Urethral Dilator Set

Used for urethral dilation.



Order Number	Reference Part Number	Dilator OD Fr	Dilator Length cm
Set			
G14185	073701	8.0–24.0	37

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

