

# Contents – Needles

## Access

Skinny Needle® with Chiba Tip . . . . .	9.2
Mitty-Pollack Needle Set. . . . .	9.3
Disposable Two-Part Trocar Needle . . . . .	9.4
Trocar Needle with Locking Stylet. . . . .	9.5
Fascial Incising Needle . . . . .	9.6


## Injection

Williams Cystoscopic Injection Needle . . . . .	9.7
---	-----



### Skinny Needle® with Chiba Tip

Used under fluoroscopy for initial location and positioning of the renal collecting system during a percutaneous procedure. The EchoTip® needle is designed to provide enhanced visualisation of the needle tip under ultrasound.



Order Number	Reference Part Number	gauge	Length cm
G14254	090010	22	20
EchoTip®			
G15290	090010-ET	22	20


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).



## Mitty-Pollack Needle Set

Used as a percutaneous entry needle to introduce wire guides. The EchoTip needle is designed to provide enhanced visualisation of the needle tip under ultrasound. The set includes an inner needle and an outer needle.



Order Number	Reference Part Number	Outer Needle gauge	Outer Needle Length cm	Inner Needle gauge	Inner Needle Length cm
<b>Set</b>					
G14266	090031	18	14	22	22
<b>EchoTip®</b>					
G15292	090031-ET	18	14	22	22


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).



## Disposable Two-Part Trocar Needle


Used for percutaneous access. The EchoTip needle is designed to provide enhanced visualisation of the needle tip under ultrasound.



Order Number	Reference Part Number	gauge	Length cm
<b>Standard</b>			
G14461	090020-S1	18	15
G14255	090020	18	20
<b>EchoTip</b>			
G15291	090020-ET	18	20


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).



## Trocar Needle with Locking Stylet


Used for introduction in nonvascular access procedures.



Order Number	Reference Part Number	gauge	Length cm
G07082	DTN-18-15.0-LS	18	15


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).



## Fascial Incising Needle

Used over a wire guide for cutting the renal fascia to facilitate dilation of a musculofascial tract.



Order Number	Reference Part Number	gauge	Length cm	Cut Width mm	Accepts Wire Guide Diameter inch
G14519	090070	18	5	4.5	0.038


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).



## Williams Cystoscopic Injection Needle

Used for cystoscopic injection into the urethra, bladder neck, and bladder wall.



Order Number	Reference Part Number	Fr	gauge	Length cm	Tip Length mm	Comments
G15296	090001-S10	3.7	23	35	8	thin-wall cannula
G16684	090001-S34	3.7	23	35	8	extra-thin-wall cannula
G14220	090001	5.0	23	35	8	
G16112	090001-S26	5.0	23	45	8	
G15276	090001-S9	5.0	25	35	8	

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).

