

Contents – Percutaneous access and drainage

Entry sets

Percutaneous Entry Set	10.2
One-Step Suprapubic Introducer	10.3

Nephrostomy catheters and sets

Cook®-Cope Loop Nephrostomy Set	10.4
Ultrathane® Nephrostomy Set with Mac-Loc® and Slip-Coat®	10.5
Universa® Loop Drainage Catheter Introductory Set	10.6
Universa® Loop Drainage Catheter Exchange Set	10.7
Percutaneous Neonatal Pigtail Nephrostomy Set	10.8
Universa® Malecot Drainage Catheter Introductory Set	10.9
Universa® Malecot Drainage Catheter Exchange Set	10.10



Suprapubic catheters and sets

Universa® Malecot Drainage Catheter - Suprapubic Set	10.11
---	-------



Percutaneous Entry Set

Used to establish a percutaneous tract into the renal pelvis for catheter placement or stone manipulation. The standard set includes a Safe-T-J® wire guide, two needles, and eight 6.0-16.0 Fr dilators.



Order Number	Reference Part Number
G14649	080000


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



One-Step Suprapubic Introducer

Intended for introducing a drainage catheter suprapubically into the bladder. The set includes a trocar stylet, a dilator, and a Peel-Away® sheath.



Percutaneous access and drainage

Order Number	Reference Part Number	Accepts Catheter* Fr	Sheath Length cm	Trocar Length cm
G25981	SANT-14-20-PLV	10.0	16.5	20
G26642	SANT-16-20-PLV	12.0	16.5	20
G26648	SANT-18-20-PLV	14.0	17.3	19.8
G26655	SANT-20-20-PLV	16.0	17	20.3
G26666	SANT-22-20-PLV	18.0	17	20.3

* The outer dimensions of catheters vary from manufacturer to manufacturer.
 Care should be taken to verify that the catheter will pass through the One-Step Peel-Away introducer prior to suprapubic insertion.
 Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.




Cook®-Cope Loop Nephrostomy Set

Used for percutaneous placement of a loop catheter in the renal pelvis for nephrostomy drainage. The standard set includes a needle, wire guides, an introducing catheter, cannulas, dilators, a retention disc, a connecting tube, and a stopcock.

CAUTION: This product contains natural rubber latex, which may cause allergic reactions in some individuals.

CAUTION: Periodic evaluation is advised; the catheter must not remain indwelling for more than 4 months.



Order Number	Reference Part Number	Fr	Length cm
Set			
G14282	085010	10.0	19.5
G14714	085012	12.0	19.5
G14917	085014	14.0	19.5
Catheter and Cannula			
G14283	085110	10.0	19.5
G14705	085112	12.0	19.5
G15214	085114	14.0	19.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.


Please see product risk information in the IFU at cookmedical.eu.



Ultrathane® Nephrostomy Set With Mac-Loc® and Slip-Coat®

Used for percutaneous placement of a loop catheter in the renal pelvis for nephrostomy drainage. The Mac-Loc provides a suture fixation system for drainage catheters. The locking cam lever provides the ability to vary the degree of loop fixation based on loop location and patient anatomy. The distal 5 cm of the catheter has a hydrophilic coating. The set includes a needle, two wire guides, a cannula, a dilator, a catheter support, a catheter, a stylet, and a connecting tube.

CAUTION: Periodic evaluation is advised; the catheter must not remain indwelling for more than 4 months.



Percutaneous access and drainage

Order Number	Reference Part Number	Fr	Length cm
G17155	085508	8.5	25
G17156	085510	10.2	25
G17157	085512	12.0	25
G17158	088514	14.0	25

Slip-Coat is a registered trademark of Argon Medical Devices BioCoatings Corp

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.


Please see product risk information in the IFU at cookmedical.eu.



Universa® Loop Drainage Catheter Introductory Set

Used to provide percutaneous urine drainage from the genitourinary system. The set includes a loop catheter, an EchoTip® needle with obturator, a curved Safe-T-J® PTFE-coated stainless steel wire guide, dilators, two straightening stylets, a connecting tube with a one-way stopcock, and a catheter support.

Percutaneous access and drainage



Order Number	Reference Part Number	Catheter Fr	Length cm	Dilators Fr	Dilator Length cm
Set					
G55802	ULCS-10-30	10.0	30	6-12	20
G55803	ULCS-12-30	12.0	30	6-14	20
G55804	ULCS-14-30	14.0	30	6-16	20
Without Stylets					
G55800	ULCS-6-15	6.0	15	5-8	20
G55801	ULCS-8-30	8.0	30	6-10	20


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Universa® Loop Drainage Catheter Exchange Set

Used to provide percutaneous urine drainage from the genitourinary system. The set includes a loop catheter, two straightening stylets, and a connecting tube with a one-way stopcock.




Percutaneous access and drainage

Order Number	Reference Part Number	Catheter Fr	Length cm
Set			
G55807	ULC-10-30	10.0	30
G55808	ULC-12-30	12.0	30
G55809	ULC-14-30	14.0	30
Without Stylets			
G55805	ULC-6-15	6.0	15
G55806	ULC-8-30	8.0	30


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Percutaneous Neonatal Pigtail Nephrostomy Set

Intended for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage in neonatal patients.



Order Number	Reference Part Number	Fr	Length cm
G14508	080106	6.0	7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.


Please see product risk information in the IFU at cookmedical.eu.



Universa® Malecot Drainage Catheter Introductory Set

Used to provide percutaneous urine drainage from the genitourinary system. The set includes a Malecot catheter, an EchoTip® needle with obturator, a curved Safe-T-J® PTFE-coated stainless steel wire guide, dilators, two straightening stylets, a connecting tube with a one-way stopcock,* and a catheter support.

**The one-way stopcock is not available on Malecot catheters that are 16 Fr and larger.*



Percutaneous access and drainage

Order Number	Reference Part Number	Catheter Fr	Length cm	Dilators Fr	Dilator Length cm
G55816	UMCS-8-30	8.0	30	6-10	20
G55817	UMCS-10-30	10.0	30	6-12	20
G55818	UMCS-12-30	12.0	30	6-14	20
G55819	UMCS-14-30	14.0	30	6-16	20
G55820	UMCS-16-30	16.0	30	6-18	20

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.




Universa® Malecot Drainage Catheter Exchange Set

Used to provide percutaneous urine drainage from the genitourinary system. The set includes a Malecot catheter, two straightening stylets, and a connecting tube with a one-way stopcock.*

**The one-way stopcock is not available on Malecot catheters that are 16 Fr and larger.*

Percutaneous access and drainage



Order Number	Reference Part Number	Catheter Fr	Length cm
G55826	UMC-8-30	8.0	30
G55827	UMC-10-30	10.0	30
G55828	UMC-12-30	12.0	30
G55829	UMC-14-30	14.0	30
G55830	UMC-16-30	16.0	30
G55831	UMC-18-30	18.0	30
G55832	UMC-20-30	20.0	30
G55833	UMC-22-30	22.0	30
G55834	UMC-24-30	24.0	30

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.


Please see product risk information in the IFU at cookmedical.eu.



Univera® Malecot Drainage Catheter Suprapubic Set

Used to provide percutaneous urine drainage from the genitourinary system. The set includes a Malecot catheter, a locking needle, and a connecting tube with a one-way stopcock.*

**The one-way stopcock is not available on Malecot catheters that are 16 Fr.*



Percutaneous access and drainage

Order Number	Reference Part Number	Catheter Fr	Length cm	Needle gage
G56334	UMCP-8-15	8.0	15	19
G55821	UMCP-8-30	8.0	30	19
G56335	UMCP-10-15	10.0	15	18
G55822	UMCP-10-30	10.0	30	18
G55823	UMCP-12-30	12.0	30	16
G55824	UMCP-14-30	14.0	30	15
G55825	UMCP-16-30	16.0	30	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

Please see product risk information in the IFU at cookmedical.eu.