

Contents – Accessories

Accessories

Foley Connector – Female	11.2
Drainage Bag Connector – Male	11.3
Tuohy-Borst Adapter	11.4


Connecting tubes

Connecting Tube	11.5
Connecting Tube – Two Female Luer Locks and One Drainage Bag Connector.	11.6
Ureteroscopy Irrigation System	11.7



Foley Connector Female

Used to attach a female Luer lock to universal taper, for connection to a Foley catheter.



Order Number	Reference Part Number
-----------------	--------------------------

G14529	053003
--------	--------


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Drainage Bag Connector Male

Used to provide connection to a drainage bag.




Order Number	Reference Part Number
G14225	050030

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.



Tuohy-Borst Adapter

Used to prevent the backflow of fluid around an instrument inserted through the working channel of a flexible or rigid ureteroscope.



Order Number	Reference Part Number	Accepts Instrumentation Fr	Comment
G18814	TBA-6	6.0	
G18813	TBAYR-6	6.0	with sidearm


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.




Connecting Tube

Used to provide connection to a drainage bag.



Order Number	Reference Part Number	Fr	Length cm
Standard			
Male Luer Lock and Drainage Bag Connector; One-Way Stopcock Attached			
G14230	050042	14.0	30
Male Luer Lock and Drainage Bag Connector; Female Luer Lock to Tuohy-Borst Adapter Attached			
G14119	050040	14.0	30


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.



Connecting Tube

Two Female Luer Locks and One Drainage Bag Connector


Used to provide connection from two catheters to a drainage bag.



Order Number	Reference Part Number	Fr	Length cm	Comments
Y-Type				
G14121	050043	14.0	30	includes two female Luer locks to Tuohy-Borst adapters


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.



Ureteroscopy Irrigation System

Used to provide controlled manual irrigation of the ureter, bladder, kidney, or uterine cavity.



Accessories

Order Number	Reference Part Number
G16918	050350

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

