

The background of the page features a central horizontal band with an abstract, marbled pattern in shades of blue, purple, and white. Above and below this band are images of medical instruments. On the left, a long, thin, silver-colored instrument with a curved tip is visible. On the right, a white, flexible, looped instrument is shown. The text 'CENTESIS & DRAINAGE PRODUCTS' is centered within the abstract band.

CENTESIS & DRAINAGE
P R O D U C T S

COOK[®]

5CM


Contents

Pneumothorax Sets	4
Wayne Pneumothorax Sets	5
Peters Pneumothorax Set	6
Landers Pneumothorax Set	7
Cook Emergency Pneumothorax Set	8
Nightengale Pneumothorax Set	9
Catheter Aspiration Set for Simple Pneumothorax	10
Simple Pneumothorax Aspiration Accessory Set	11
Thal-Quick Chest Tube Sets and Trays	12-13
Thal-Quick Chest Tube Adapter	14
Thoracentesis Set	15
Pericardiocentesis Set	16
Lock Pericardiocentesis Set	17
Tilkian Pericardiocentesis Tray	18
Teitel Pediatric Pericardiocentesis Set	19
Fuhrman Pleural/Pneumopericardial Drainage Set	20
Pleural/Pneumopericardial Drainage Set	21
Emergency Pericardiocentesis Set and Tray	22
Peritoneal Lavage Set and Tray	23
References	24
Acknowledgements	24
Videotapes	24

PNEUMOTHORAX SETS



Used to treat spontaneous and tension pneumothorax. The adapter on the tubing can be connected to a suction apparatus or the Cook Chest Drain Valve. Supplied sterile in peel-open packages. Intended for one-time use.




SET ORDER NUMBER	CATHETER			NEEDLE		CONNECTING TUBE			Global Product Number
	French Size	Length	Side- ports	Gage	Length	Point	Length	Remarks	
Pediatric									
C-TPTP-100	6.3	19 cm	6	21	12 cm	trocar	20 cm	Catheter is radiopaque TFE.	G03465
Standard									
C-TPT-100	9.0	29 cm	20	18	20 cm	trocar	30 cm	Catheter is radiopaque TFE.	G03301
C-TPT-200	10.2	29 cm	20	18	20 cm	trocar	30 cm	Catheter is radiopaque polyurethane.	G03302

WAYNE PNEUMOTHORAX SETS



Used for the relief of simple, spontaneous, iatrogenic and tension pneumothorax. Supplied sterile in peel-open packages. Intended for one-time use.




SET ORDER NUMBER	CATHETER			NEEDLE OBTURATOR			CONNECTING TUBE	Global Product Number
	French Size	Length	Side- ports	Gage	Length	Point	Length	
C-UTPT-1020-WAYNE	10.2	29 cm	20	19	16.5 cm	trocar	30 cm	G09922
C-UTPT-1400-WAYNE	14.0	29 cm	20	19	16.5 cm	trocar	30 cm	G07910

PETERS PNEUMOTHORAX SET



Used to treat spontaneous and tension pneumothorax. Larger catheter sideports also permit effective fluid drainage. The adapter on the tubing can be connected to a suction apparatus or the Cook Chest Drain Valve. Supplied sterile in peel-open packages. Intended for one-time use.



SET ORDER NUMBER	CATHETER			NEEDLE			CONNECTING TUBE	Global Product Number
	French Size	Length	Side- ports	Gage	Length	Point	Length	
C-TPT-100-PTR	9.0	20 cm	20	18	20 cm	trocár	15 cm	G05299


LANDERS PNEUMOTHORAX SET



Used to treat spontaneous and tension pneumothorax.

- The adapter on the tubing can be connected to a suction apparatus or the Cook Chest Drain Valve.
- Introduction is performed by blunt dissection with use of blunt obturator cannula.

Supplied sterile in peel-open packages. Intended for one-time use.



SET ORDER NUMBER	CATHETER			BLUNT TIP OBTURATOR		CONNECTING TUBE	Global Product Number
	French Size	Length	Side- ports	Gage	Length	Length	
C-TPT-100-LDR	9.0	20 cm	19	18	23.4 cm	30 cm	G05798


COOK EMERGENCY PNEUMOTHORAX SET




Used to treat simple and tension pneumothorax. Supplied sterile in peel-open packages. Intended for one-time use.




CATHETER INTRODUCER NEEDLE
Reinforced radiopaque FEP catheter
and stainless steel, appropriate length needle




COOK CHEST DRAIN VALVE
Latex free




CONNECTING TUBE
Polyvinylchloride



ONE-WAY
STOPCOCK
Plastic



SYRINGE



MOLNAR DISC
WITH PULL TIE

SET ORDER NUMBER ¹	CATHETER		NEEDLE		CONNECTING TUBE	Global Product Number
	French Size	Length	Gage	Point	Length	
C-TPTS-8.5-6.0-FSNS	8.5	6.0 cm	15	lancet	30 cm	G13722


¹Set consists of items shown above and the following items: alcohol prep, povidone-iodine swabstick and porous transparent tape.

NIGHTENGAL PNEUMOTHORAX SET



Used to treat spontaneous and tension pneumothorax in emergency situations.

- Set is preassembled to expedite access procedure and packaged appropriately for emergency transport.
 - Three-way stopcock permits manual evacuation of air without disconnecting from Cook Chest Drain Valve if indicated.
- Supplied sterile in peel-open packages. Intended for one-time use.




SET ORDER NUMBER ¹	CATHETER			NEEDLE			Global Product Number
	French Size	Length	Side- ports	Gage	Length	Point	
C-TPT-100-NTG	9.0	21.5 cm	20	18	20 cm	trocar	G05282

¹Set consists of items shown above and the following items: 25 gage needle, lidocaine, povidone-iodine swabsticks, disposable syringe and disposable scalpel.

CATHETER ASPIRATION SET FOR SIMPLE PNEUMOTHORAX



Used to treat simple pneumothorax. Use of standard percutaneous entry (Seldinger) technique facilitates controlled, atraumatic catheter introduction. Supplied sterile in peel-open packages. Intended for one-time use.




SET ORDER NUMBER	CATHETER			WIRE GUIDE		NEEDLE		Global Product Number
	French Size	Length	Side- ports	Diameter	Length	Gage	Length	
C-CASP-901-FORD	9.0	25 cm	10	.038 inch	50 cm	18	7 cm	G05296

SIMPLE PNEUMOTHORAX ASPIRATION ACCESSORY SET



Used in conjunction with pneumothorax aspiration catheters if incomplete expansion of the lung persists after mechanical aspiration of air. Connection of the one-way Cook Chest Drain Valve to the catheter permits continued re-expansion of the lung by natural forces, preventing inadvertent backflow of air into the pleural space. Supplied sterile in peel-open packages. Intended for one-time use.


**SET
ORDER
NUMBER**
Remarks
**Global
Product
Number**
C-CASP-A-FORD

SET consists of items shown above.

G05297


THAL-QUICK CHEST TUBE SETS AND TRAYS



Used for pleural fluid drainage. Chest tube placement is facilitated utilizing standard percutaneous entry (Seldinger) technique.

- The chest tube inserter, fitted to the chest tube, provides a smooth transition for ease of insertion over the wire guide.
- Sideports are positioned at the distal end of the chest tube.

Supplied sterile in peel-open packages. Intended for one-time use.




THAL-QUICK CHEST TUBE SETS AND TRAYS

ORDER NUMBER	THAL-QUICK CHEST TUBE			Global Product Number
	French Size	Length	Sideports	
Sets ¹				
C-TQTS-800	8.0	18 cm	3	G06786
C-TQTS-1000	10.0	20 cm	3	G06787
C-TQTS-1200	12.0	22 cm	3	G05161
C-TQTS-1400	14.0	41 cm	4	G05777
C-TQTS-1600	16.0	41 cm	4	G04891
C-TQTS-1800	18.0	41 cm	4	G05453
C-TQTS-2000	20.0	41 cm	4	G05454
C-TQTS-2400	24.0	41 cm	4	G04137
C-TQTS-2600	26.0	41 cm	4	G04138
C-TQTS-2800	28.0	41 cm	4	G04139
C-TQTS-3200	32.0	41 cm	4	G07243
C-TQTS-3600	36.0	41 cm	4	G07241
Trays ²				
C-TQTSY-800	8.0	18 cm	3	G06791
C-TQTSY-1000	10.0	20 cm	3	G06885
C-TQTSY-1200	12.0	22 cm	3	G05464
C-TQTSY-1400	14.0	41 cm	4	G07090
C-TQTSY-1600	16.0	41 cm	4	G05462
C-TQTSY-1800	18.0	41 cm	4	G05461
C-TQTSY-2000	20.0	41 cm	4	G05463
C-TQTSY-2400	24.0	41 cm	4	G04220
C-TQTSY-2600	26.0	41 cm	4	G04225
C-TQTSY-2800	28.0	41 cm	4	G04226
C-TQTSY-3200	32.0	41 cm	4	G07242
C-TQTSY-3600	36.0	41 cm	4	G07244


¹ All sets include the items shown on the previous page.

² All trays consist of items shown on the previous page and the following items: 25 gage needle, 22 gage needle, lidocaine, povidone-iodine swabsticks, fenestrated drape, gauze sponges, disposable syringes, suture with needle, povidone-iodine ointment, CSR wrap, prep tray and needle holder cup.

THAL-QUICK CHEST TUBE ADAPTER




Used to provide separate sampling or infusion port within Thal-Quick Chest Tubes. Supplied sterile in peel-open packages. Intended for one-time use.




ORDER NUMBER	Accepts	Global Product Number
C-TQF-28.0	24.0 thru 36.0 French	G07467

THORACENTESIS SET



Used to remove fluid from the pleural space. Supplied sterile in peel-open packages. Intended for one-time use.




SET ORDER NUMBER	CATHETER			STIFFENING CANNULA		WIRE GUIDE			NEEDLE		Global Product Number
	French Size	Length	Side- ports	Gage	Length	Diameter	Length	Tip	Gage	Length	
C-THS-850	8.5	22 cm	4	18	24 cm	.038 inch	70 cm	3 mm J	18	10 cm	G03286

PERICARDIOCENTESIS SET




Used to drain fluid from the pericardial sac. The “alligator clip” cable is attached to the hub of the needle for constant electrocardiographic monitoring during needle insertion. Supplied sterile in peel-open packages. Intended for one-time use.




SET ORDER NUMBER	CATHETER			STIFFENING CANNULA		WIRE GUIDE			NEEDLE		Global Product Number
	French Size	Length	Side- ports	Gage	Length	Diameter	Length	Tip	Gage	Length	
C-PCS-850	8.5	22 cm	4	18	24 cm	.038 inch	70 cm	3 mm J	18	15 cm	G03282

LOCK PERICARDIOCENTESIS SET



Used to drain fluid from the pericardial sac. Supplied sterile in peel-open packages. Intended for one-time use.




SET ORDER NUMBER	CATHETER			WIRE GUIDE			NEEDLE		Global Product Number
	French Size	Length	Side- ports	Diameter	Length	Tip	Gage	Length	
C-PCS-830-LOCK	8.3	40 cm	6	.038 inch	70 cm	3 mm J	18	7 cm	G02433

TILKIAN PERICARDIOCENTESIS TRAY




Used to drain fluid from the pericardial sac. The “alligator clip” cable is attached to the hub of the needle for constant electrocardiographic monitoring during needle insertion. Supplied sterile in peel-open packages. Intended for one-time use.




TRAY ORDER NUMBER ¹	CATHETER			WIRE GUIDE			NEEDLE		Global Product Number
	French Size	Length	Side- ports	Diameter	Length	Tip	Gage	Length	
C-PCSY-830-AT	8.3	40 cm	6	.038 inch	70 cm	straight and 3 mm J	18	8 cm	G05252

¹ Tray consists of items shown above and the following items: 25 gage needle, 22 gage needle, lidocaine, povidone-iodine swabsticks, fenestrated drape, gauze sponges, disposable syringes, disposable scalpel, suture with needle, povidone-iodine ointment, CSR wrap, prep tray and needle holder cup.

TEITEL PEDIATRIC PERICARDIOCENTESIS SET



Used to drain fluid from the pericardial sac. The “alligator clip” cable is attached to the hub of the needle for constant electrocardiographic monitoring during needle insertion. Supplied sterile in peel-open packages. Intended for one-time use.



SET ORDER NUMBER	CATHETER			WIRE GUIDE			NEEDLE		Global Product Number
	French Size	Length	Side- ports	Diameter	Length	Tip	Gage	Length	
C-PCS-500-TTL	5.0	30 cm	6	.028 inch	60 cm	3 mm J	19	6 cm	G05251

FUHRMAN PLEURAL/PNEUMOPERICARDIAL DRAINAGE SET



Used for evacuation of air from the pericardial sac or to drain air or fluid from the pleural space. Supplied sterile in peel-open packages. Intended for one-time use.



INTRODUCER NEEDLE
Stainless steel



WIRE GUIDE
TFE-coated stainless steel wire guide



DILATOR
8.0 French radiopaque 12 cm long



PIGTAIL CATHETER
Radiopaque polyurethane




MULTIPURPOSE TUBING ADAPTER
Plastic

SET ORDER NUMBER	CATHETER			WIRE GUIDE			NEEDLE		Global Product Number
	French Size	Length	Side- ports	Diameter	Length	Tip	Gage	Length	
C-PPD-850	8.5	15 cm	6	.028 inch	40 cm	2 mm J	18	7 cm	G03974

PLEURAL/PNEUMOPERICARDIAL DRAINAGE SET



Used for evacuation of fluid from the pleural space in neonates.
Supplied sterile in peel-open packages. Intended for one-time use.



SET ORDER NUMBER ¹	CATHETER			NEEDLE		Global Product Number
	French Size	Length	Side- ports	Gage	Length	
C-PPD-500-WOOD	5.0	8 cm	10	18	16 cm	G11889
C-PPD-600-WOOD	6.0	8 cm	10	18	16 cm	G11890

¹Set consists of items shown above and the following items: Disposable scalpel, needle holder, fenestrated drape and betadine solution.


EMERGENCY PERICARDIOCENTESIS SET AND TRAY



Used to drain fluid from the pericardial sac.

- Coaxial catheter needle allows one-step emergency access.
- The “alligator clip” cable is attached to the hub of the needle for constant electrocardiographic monitoring during insertion.

Supplied sterile in peel-open packages. Intended for one-time use.




ORDER NUMBER	CATHETER				NEEDLE		CONNECTING TUBE	Global Product Number
	French Size	Length	Side- ports	Sideport Diameter	Gage	Length	Length	
Set								
C-PCS-1000-SPAC	10.0	19 cm	2	.065 inch	18	21 cm	7.5 cm	G05250
Tray ¹								
C-PCSY-1000-SPAC	10.0	19 cm	2	.065 inch	18	21 cm	7.5 cm	G05298

¹ Tray consists of items shown above and the following items: 25 gage needle, 22 gage needle, lidocaine, povidone-iodine swabsticks, fenestrated drape, gauze sponges, disposable syringe, suture with needle, povidone-iodine ointment, CSR wrap, prep tray and needle holder cup.

PERITONEAL LAVAGE SET AND TRAY



Used to infuse and withdraw rinsing solutions in the peritoneal cavity. Supplied sterile in peel-open packages. Intended for one-time use.



ORDER NUMBER	CATHETER			WIRE GUIDE		NEEDLE		Global Product Number
	French Size	Length	Side- ports	Diameter	Length	Gage	Length	
	Set							
C-PLS-901	9.0	20 cm	90	.038 inch	40 cm	18	7 cm	G05267
Tray ¹								
C-PLSY-901	9.0	20 cm	90	.038 inch	40 cm	18	7 cm	G05268

¹ Tray consists of items shown above and the following items: 25 gage needle, 22 gage needle, lidocaine, povidone-iodine swabsticks, fenestrated drape, gauze sponges, disposable syringes, disposable scalpel, suture with needle, povidone-iodine ointment, CSR wrap, prep tray and needle holder cup.

REFERENCES

- Sargent, E. N., Turner, A. F.: "Emergency Treatment of Pneumothorax: A Simple Catheter Technique for Use in the Radiology Department," *American Journal of Roentgenology*, 109 (1970), 531-535.
- Lock, J. E., Bass, J. L., Kulik, T. J., Fuhrman, B. P.: "Chronic Percutaneous Pericardial Drainage with Modified Pigtail Catheters in Children," *American Journal of Cardiology*, 53 (1984) 1179-1182.
- Fuhrman, B. P., Landrum, B. G., Ferrara, T. B., Steinhorn, D. M., Connell, A. P., Smith-Wright, D. L., Green, T. P.: "Pleural Drainage Using Modified Pigtail Catheters," *Critical Care Medicine*, 14 (1986), 575-576.
- Tilkian, A., Daily, E. K.: "Pericardiocentesis and Drainage," Chapter 10, pgs. 231-256 in *Cardiovascular Procedures: Diagnostic Techniques and Therapeutic Procedures*. The C. V. Mosby Company, St. Louis, Missouri, 1986.
- Thal, A. P., Quick, K. L.: "A Guided Chest Tube for Safe Thoracostomy," *Surgery, Gynecology and Obstetrics*, 167 (1988), 517.
- Vallee, P., Sullivan, M., Richardson, H., Bivins, B., Tomlanovitch, M.: "Sequential Treatment of a Simple Pneumothorax," *Annals of Emergency Medicine*, 17 (1988), 936-942.
- Conces, D. J., Jr., Tarver, R. D., Gray, W. C., Pearcy, E. A.: "Treatment of Pneumothoraces Utilizing Small Caliber Chest Tubes," *Chest*, 94 (1988), 55-57.
- Wood, B., Dubik, M.: "A New Device for Pleural Drainage in Newborn Infants," *Pediatrics*, 96 (1995), 955-956.
- Ahmed, M. Y., Silver, P., Nimkoff, L., Sagy, M.: "The Needle-Wire-Dilator Technique for the Insertions of Chest Tubes in Pediatric Patients," *Pediatric Emergency Care*, 11 (1995), 252-254.
- Higginson, I.: "Emergency Pneumothorax Set," *Pre-Hospital Immediate Care*, 2 (1998), 104-105.
- Gammie, J. S., Banks, M. C., Fuhrman, C. R., Pham, S. M., Griffith, B. P., Keenan, R. J., Luketich, J. D.: "The Pigtail Catheter for Pleural Drainage: A Less Invasive Alternative to Tube Thoracostomy," *Journal of the Society of Laparoendoscopic Surgeons*, 3 (1999), 57-61.

ACKNOWLEDGEMENTS

- Landers, C. F., M.D., Division of Pulmonary and Critical Care Medicine, D. N. Sharp Memorial Hospital, San Diego, California.
- Teitel, D., M.D., Department of Pediatric Cardiology, University of California at San Francisco, San Francisco, California.

VIDEOTAPES

VHS (1/2 inch) videotapes are available. Allow 1-2 weeks for delivery. For listing of updates or additional information on available videotapes, contact COOK CRITICAL CARE.

ORDER NUMBER

C-VC-16 PERCUTANEOUS INSERTION OF THE THAL-QUICK CHEST TUBE

AUTHORS: Thal, A. P., M.D., F.A.C.S., Quick, K. L.,
B.S.N., Department of Surgery, Northeastern
Regional Hospital, Las Vegas, New Mexico.
(Duration: 4 min.)

ORDER NUMBER

C-VC-39 COOK PNEUMOTHORAX CATHETER SET WAYNE PNEUMOTHORAX CATHETER SET

AUTHOR: Wayne, M. A., M.D., F.A.C.E.P., Medical
Director, Whatcom County Emergency Medical
Service, Bellingham, Washington.
(Duration: 5 min.)

COOK®

Interventional Critical Care Products

www.cookcriticalcare.com

Distributed by:

COOK MEDICAL INCORPORATED
P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.
Phone: 812 339-2235, Toll Free: 800 457-4500, Toll Free Fax: 800 554-8335

COOK (CANADA) INC.
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905 640-7110, Toll Free: 800 668-0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.
Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86