

Greene Renal Transplant Stent Set

- [Catalog Page](#)
- [Component List](#)
- [Resource Library](#)



Used to establish temporary internal drainage from the ureteropelvic junction to the bladder following renal transplant. Non-radiopaque ink bands serve as a visual indicator so the physician can anticipate where the coil(s) will form as the stent is placed. Set Includes: Stent, Wire Guide, Catheter, and Stent Positioner

CAUTION: Periodic evaluation is advised; stent must not remain indwelling more than 12 months.

VZP	Objednací číslo		Fr	Délka stentu
0052616	U-435520	NS	5.0	8 to 20 cm
0052616	U-435620	NS	6.0	8 to 20 cm
	U-435632		6.0	22 to 32 cm
0052616	U-435720	NS	7.0	8 to 20 cm
0052616	U-435820	NS	8.5	8 to 20 cm
	U-435832		8.5	22 to 32 cm

Michaela Gray

Mobil: 602 561 576

strana: 1 / 1



Sídlo:

Na Dolinách 128/36
147 00 Praha 4
tel.: 241 431 122
fax: 241 430 019

Provozovna:

Vlastibořská 2790/4
193 00 Praha 9
tel.: 296 328 328
fax: 296 328 320

IČ: 14892901

DIČ: CZ14892901

Bankovní spojení:

ČSOB Praha

č.ú.: 758073/0300

Registrace v obchodním

rejstříku Městského soudu

Praha

Odd. C, vložka 1299

www.medial.cz