




KOYLE DIAPER STENTS

Used for stenting the urethra after hypospadias or epispadias repair and to provide postoperative drainage of the bladder. The externalized portion of the stent is of sufficient length to allow drainage into a second or outside diaper, keeping the inside diaper dry during the healing process. Supplied sterile in peel-open packages. Intended for one-time use.



ORDER NUMBER	French Size	Length
STENT ONLY		
KDSCF-060050	6.0	50 cm
KDSCF-080050	8.0	50 cm
STENT AND INTRODUCER		
KDSCFS-060050	6.0	50 cm
KDSCFS-080050	8.0	50 cm

Mgr. Michaela Gray
Mob: 602 561 576