




Williams

CYSTOSCOPIC INJECTION NEEDLE




Inject with confidence.



Williams

CYSTOSCOPIC INJECTION NEEDLE



The Williams Cystoscopic Injection Needle is used for the injection of a clinically approved therapeutic agent for indications in the urethra, bladder neck and bladder wall.

Order Number	Reference Part Number	Fr	Needle gage	Length cm	Tip Length mm
G16684	090001-S34	3.7	23	35	8
G14220	090001	5.0	23	35	8
G16112	090001-S26	5.0	23	45	8

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
E-mail: cau.custserv@cookmedical.com

