



Endoureterotomy Stent Set

Ultrathane®

Used for temporary internal drainage from the ureteropelvic junction to the bladder following incision of a stricture. The segment with no sideports prevents postoperative narrowing of the ureteral lumen while preventing ingrowth of the ureteral wall into the stent. Set Includes: Stent, Release Sleeve, and Stent Positioner

CAUTION: Periodic evaluation is advised; stent must not remain indwelling more than eight (8) weeks.

Global Product Number	Order Number	French Size	Proximal Segment Length	Overall Length
G17144	UESS-06951318	6.0 to 9.5	13 cm	18 cm
G17143	UESS-06951320	6.0 to 9.5	13 cm	20 cm
G17142	UESS-06951322	6.0 to 9.5	13 cm	22 cm
G17141	UESS-06951324	6.0 to 9.5	13 cm	24 cm
G17140	UESS-06951326	6.0 to 9.5	13 cm	26 cm
G17139	UESS-06951328	6.0 to 9.5	13 cm	28 cm
G17138	UESS-06951330	6.0 to 9.5	13 cm	30 cm
G17136	UESS-07141320	7.0 to 14.0	13 cm	20 cm
G17135	UESS-07141322	7.0 to 14.0	13 cm	22 cm
G17133	UESS-07141324	7.0 to 14.0	13 cm	24 cm
G17132	UESS-07141326	7.0 to 14.0	13 cm	26 cm
G17131	UESS-07141328	7.0 to 14.0	13 cm	28 cm
G17130	UESS-07141330	7.0 to 14.0	13 cm	30 cm



Siedmon Antegrade AQ® Stent Set

Used for temporary internal drainage from the ureteropelvic junction to the bladder. Set provides for antegrade placement of a ureteral double pigtail stent. The AQ® coating is a microthin layer of hydrophilic polymer that when activated, attracts and holds water and other liquids to the wire guide, creating a low-friction surface. Set includes: stent and introducer

CAUTION: Periodic evaluation is advised; stent must not remain indwelling for more than six (6) months.

Global Product Number	Order Number	French Size	Stent Length
G19119	025306-S4-AQ	6.0	50 cm