

## C-Flex® Multi-Length Ureteral Stent Set



Used for temporary internal drainage from the ureteropelvic junction to the bladder. Set includes stent, wire guide and catheter.

**CAUTION: Periodic evaluation is advised; stent must not remain indwelling more than 6 months.**

C-Flex® Multi-Length Ureteral Stent Set

VZP 0051475

Objednací kód	Fr	Délka stentu cm	Komentář	
U-036300-T-8-20	3.7	8 to 20	with tether	S
U-036300	3.7	21 to 32 cm		NS

  

INCLUDED COMPONENTS	
Description	Quantity
Multi-Length Stent – Radiopaque	1
Straight Safety Wire Guide .025 inch (0.64 mm) Stainless Steel 145 cm long with 3 cm flexible tip	1
Catheter 4.0 French Radiopaque polyurethane 70 cm long	1
Stent Positioner Radiolucent 3.7 French 50 cm long	1

strana: 1 / 1



**Sídlo:**  
Na Dolinách 128/36  
147 00 Praha 4  
tel.: 241 431 122  
fax: 241 430 019

**Provozovna:**  
Vlastibořská 2790/4  
193 00 Praha 9  
tel.: 296 328 328  
fax: 296 328 320

IČ: 14892901  
DIČ: CZ14892901  
Bankovní spojení:  
ČSOB Praha  
č.ú.: 758073/0300

Registrace v obchodním  
rejstříku Městského soudu  
Praha  
Odd. C, vložka 1299  
www.medial.cz